CooperSurgical®

## LEEP PRECISION™ Integrated System

120 VAC Model LP-10-120

230 VAC Model LP-10-220

**Operating Manual (English)**

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---

# LEEP PRECISION™ Integrated System, 120 VAC and 230 VAC

Models LP-10-120 and LP-10-220

## Table of Contents

Section	Content	Page
<b>1.</b>	<b>Description</b> .....	<b>5</b>
1.1	Introduction .....	5
1.2	LEEP PRECISION Generator Description .....	5
1.3	LEEP PRECISION Smoke Evacuator Description .....	5
1.4	LEEP PRECISION Cart (with attached LEEP PRECISION Integration Unit) Description .....	5
<b>2.</b>	<b>Unpacking and Assembly</b> .....	<b>6</b>
2.1	Unpacking Shipping Carton .....	6
2.2	Unpacking and Installing the LEEP PRECISION Smoke Evacuator .....	7
2.3	Unpacking and Installing the LEEP PRECISION Generator .....	7
2.4	Installing the Filters and Tubing to the LEEP PRECISION Smoke Evacuator .....	8
2.5	Connecting the LEEP PRECISION Integrated System to the Wall Outlet .....	8
2.6	Installing the Foot Pedal Switch on the LEEP PRECISION Generator .....	9
2.7	Installing the Active and Dispersive Electrodes .....	9
<b>3.</b>	<b>LEEP PRECISION Generator Features</b> .....	<b>10</b>
<b>4.</b>	<b>Front and Rear Panels of the LEEP PRECISION Generator and Integration Unit</b> .....	<b>12</b>
4.1	Front Panel of the LEEP PRECISION Generator .....	12
4.2	Rear Panel of the LEEP PRECISION Generator .....	13
4.3	LEEP PRECISION Integration Unit .....	13
<b>5.</b>	<b>Professional Use Guide</b> .....	<b>14</b>
5.1	LEEP PRECISION Cart Use .....	14
5.2	Indications .....	14
5.3	Contraindications .....	15
5.4	LEEP Procedure and Technique .....	15
5.5	Safety Precautions .....	15
5.6	Electrosurgical Procedures .....	16
<b>6.</b>	<b>Electrosurgical Precautions</b> .....	<b>19</b>
<b>7.</b>	<b>LEEP PRECISION Patient Return Electrode</b> .....	<b>21</b>
<b>8.</b>	<b>Turning on the LEEP PRECISION Integrated System</b> .....	<b>22</b>
8.1	Practice .....	22
8.2	The Power Setting .....	24
8.3	Cutting Techniques .....	24
8.4	Criteria of a Good Cutting Technique .....	25
8.5	Coagulating .....	25
8.6	Coagulating Technique .....	25
8.7	Technique Guidance .....	26
<b>9.</b>	<b>Powering Down the LEEP PRECISION Integrated System</b> .....	<b>26</b>
<b>10.</b>	<b>Maintenance</b> .....	<b>27</b>
<b>11.</b>	<b>Accessories</b> .....	<b>27</b>
11.1	Disposable Hand Switch Pencils .....	27

## Table of Contents (continued)

Section	Content	Page
<b>12.</b>	<b>Troubleshooting</b> .....	<b>28</b>
<b>13.</b>	<b>Specifications</b> .....	<b>29</b>
13.1	LEEP PRECISION Integrated System .....	29
13.2	LEEP PRECISION Generator .....	29
13.3	LEEP PRECISION Smoke Evacuator .....	30
13.4	LEEP PRECISION Cart with Integration Unit .....	30
13.5	Power Output Characteristics .....	31
13.6	Power Output at Various Load Resistance .....	32
<b>14.</b>	<b>LEEP PRECISION Smoke Evacuator</b> .....	<b>33</b>
14.1	System Description .....	33
14.2	Initial System Setup .....	33
14.3	ULPA Filter Installation .....	33
14.4	Operating Instructions .....	33
14.5	Maintenance Procedures .....	34
14.6	Cleaning and Disinfection .....	34
14.7	Cautions .....	34
<b>15.</b>	<b>EMC Compliance Information for the LEEP PRECISION Integrated System</b> .....	<b>35</b>
<b>16.</b>	<b>Guidance and Manufacturer's Declaration – Electromagnetic Immunity</b> .....	<b>36</b>
<b>17.</b>	<b>Liability Statement</b> .....	<b>38</b>
<b>18.</b>	<b>Warranty</b> .....	<b>39</b>
<b>19.</b>	<b>Service and Repair</b> .....	<b>39</b>
<b>20.</b>	<b>Explanation of Symbols</b> .....	<b>40</b>

CooperSurgical®

95 Corporate Drive  
Trumbull, CT 06611 USA  
Phone: (800) 243-2974  
Fax: (800) 262-0105  
[www.coopersurgical.com](http://www.coopersurgical.com)  
International

Phone: +1 (203) 601-9818  
Fax: +1 (203) 601-4747

CooperSurgical, Inc.  
95 Corporate Drive  
Trumbull, CT 06611 USA

# Section 1 Description

## 1.1 Introduction

The LEEP PRECISION™ Integrated System, which includes the following components:

LEEP PRECISION Integrated System				
LP-10-120	Generator, 120 VAC	Smoke Evacuator 120 VAC	Cart	Integration Unit, 120 VAC
LP-10-220	Generator, 230 VAC	Smoke Evacuator 230 VAC	Cart	Integration Unit, 230 VAC

## 1.2 LEEP PRECISION Generator Description

The **LEEP PRECISION Generator** has the following features:

- Isolated power output and LED display located in the front for power selection, delivery, and ease of use
- Flush faceplate membrane to facilitate operation and cleaning
- Microprocessor-based controls for increased precision, accuracy, reproducibility, and safety
- A choice of CUT, BLEND, and COAG waveforms to accommodate subtle differences in technique and electrode performance
- Hand or foot activation for precise operation
- Patient return monitor that disables power automatically if Patient Return Electrode is loose
- Self-monitors all detectors

## 1.3 LEEP PRECISION Smoke Evacuator Description

The CooperSurgical **LEEP PRECISION Smoke Evacuator System** three-stage air filtration system is used to remove airborne particulate plumes produced during office and surgical procedures, and has the following features:

- Low noise level
- Triple filtration of air provides an efficiency level of 0.014 microns rated at 99.999 percent. This includes a pre-filter, a charcoal filter for odor removal, and a final safety filter placed after the charcoal filter
- High air flow for plume collection
- Attaches to the CooperSurgical LEEP PRECISION Integrated System

## 1.4 LEEP PRECISION Cart (with attached LEEP PRECISION Integration Unit) Description

**The LEEP PRECISION Cart:**

- The LEEP PRECISION Generator and Smoke Evacuator in one unit
- Heavy-duty casters for mobility
- Storage shelves

**Purpose of the LEEP PRECISION Integration Unit:**

- The LEEP PRECISION Integration Unit controls the ON/OFF power of the entire LEEP PRECISION Integrated System. (The LEEP PRECISION Generator and Smoke Evacuator are powered through the LEEP PRECISION Integration Unit.)

## Section 2 Unpacking and Assembly

- Locate all product boxes.

### 2.1 Unpacking Shipping Carton

1. Remove the LEEP PRECISION Smoke Evacuator.
2. Remove the LEEP PRECISION Generator box.
3. Remove the LEEP PRECISION Cart.

**CAUTION: DO NOT PICK THE CART UP BY ITS HANDLES TO REMOVE FROM THE SHIPPING BOX. THE HANDLES ARE NOT DESIGNED TO ACCOMMODATE THE FULL WEIGHT OF THE CART.**


Photo A


Figure 1

## 2.2 Unpacking and Installing the LEEP PRECISION Smoke Evacuator

After unpacking the LEEP PRECISION Smoke Evacuator:


1. Undo the Velcro® strap located at the rear of the LEEP PRECISION Cart.
2. Plug in the connector [located at (B) in Figure 1] from the LEEP PRECISION Integration Unit to the rear of the LEEP PRECISION Smoke Evacuator.
3. Put the LEEP PRECISION Smoke Evacuator in the designated area (slide in place).
4. Refasten the Velcro strap.

Install the LEEP PRECISION Generator next.

## 2.3 Unpacking and Installing the LEEP PRECISION Generator

After unpacking the LEEP PRECISION Generator:

1. Turn on the power switch on the LEEP PRECISION Generator (put the switch in the ON position ("I")).  
The LEEP PRECISION Generator can now be placed into the sleeve which hangs over the LEEP PRECISION Smoke Evacuator.
2. Rotate the cart so the handle is in the front.  
The plug on the LEEP PRECISION Integration Unit (inside the sleeve) fits into the receptacle of the LEEP PRECISION Generator. Gently press the LEEP PRECISION Generator into place.

## 2.4 Installing the Filters and Tubing to the LEEP PRECISION Smoke Evacuator

### 2.4.1 ULPA Filter Installation

Tilt the LEEP PRECISION Smoke Evacuator forward and insert the large ULPA Filter Cylinder with the Air Flow arrow pointing downward (see Photo B).

Photo B

Photo C

### 2.4.2 Installing the Pre-Filter

Prepare the unit for use by inserting a clean, disposable Pre-Filter (Single Patient Use) onto the ULPA Filter Cylinder (see Photo C). Be sure the device is firmly seated.

### 2.4.3 Connecting the Tubing

There are two tubing hookup options:

1. For procedures requiring close-proximity plume removal (i.e., vaginal speculum)
2. For procedures requiring open-area plume removal (i.e., external lesions)

#### For procedures requiring close-proximity plume removal (i.e., vaginal speculum)

Assemble the  $\frac{3}{8}$ -inch Reducer (REF 6083) onto the port on the disposable Pre-Filter (REF 6081) top with a slight twisting motion. Attach one end of an appropriate length of  $\frac{3}{8}$ -inch ID Evacuation Tubing (REF 6084) to the Reducer connector and direct the other end to the patient and any appropriate device being used, such as a vaginal speculum equipped with a Smoke Evacuation Adapter (see Photo D).


Photo D

#### For procedures requiring open-area plume removal (i.e., external lesions)

Assemble the sterile Disposable Evacuation Tubing (REF 6084) directly into the top of the Pre-Filter. Position the opposite end over the site to be treated.

## 2.5 Connecting the LEEP PRECISION Integrated System to the Wall Outlet

Make sure the ON/OFF Switch **12** for the LEEP PRECISION Integrated System, located on the LEEP PRECISION Integration Unit, is in the OFF position ("O"). Next, attach the power cord to the receptacle **13** on the LEEP PRECISION Integration Unit and then into a hospital-grade grounded wall outlet to obtain power [refer to (A) in Figure 1].

## 2.6 Installing the Foot Pedal Switch on the LEEP PRECISION Generator

Connect the Foot Pedal Switch to the socket **5** shown in Figure 2 if foot control is desired.


Figure 2

## 2.7 Installing the Active and Dispersive Electrodes

### 2.7.1 Placement of the Active Electrode

- Connect the Active Electrode to socket **7** shown in Figure 2.

#### 2.7.1.1 For Foot Switch-Operated Electrode Use (with Adapter)

A 4 mm Reusable Electrosurgical Adapter is supplied with the LEEP PRECISION Integrated System (found in the LEEP PRECISION Generator box). Use the Foot Switch-Operated Electrode with the Adapter. The Electrode is plugged into the Adapter, which is then plugged into socket **7** shown in Figure 2.

### 2.7.2 Placement of the Dispersive Electrode or Patient Return Electrode

When using an electrosurgical system, it is very important that all of the current delivered to the patient returns correctly to the LEEP PRECISION Generator via the Patient Return Electrode only.

- Connect the Dispersive Electrode to socket **8**. Refer to Figure 2.
- The patient must be positioned correctly on the operating table. The patient and operator must not come in contact with any metal conductive surfaces.
- The Patient Return Electrode must securely contact a vascular area close to the operating site. For a gynecology procedure, the preferred site is the patient's thigh. The contact area must be clean, free of body lotions, shaved, and massaged for good circulation. The Patient Return Electrode contact area must be maximized and frequently checked for uniform contact during the procedure, especially if the patient has moved or if liquids have come in contact with the Patient Return Electrode. The Patient Return Electrode must NEVER be placed so as to allow the patient's heart to be in the pathway from the active electrode.
- Power delivery to the operative site may be decreased appreciably if alternate pathways exist; for example, through the metal operating table, crossed Active Electrode/Patient Return Electrode Cables, etc.

Figures 3 through 5 show the proper and improper ways of hooking up and using the various Active Electrodes and Patient Return Electrodes on the patient.

**PROPER**


Figure 3

**IMPROPER**


Figure 4

**IMPROPER**


Figure 5

### 2.7.2.1 Patient Return Electrode Use

There are two varieties of Patient Return Electrodes that can be used with the LEEP PRECISION Integrated System:

1	2
<b>Patient Return Electrode With Connector</b>	<b>Patient Return Electrode Without Connector</b>
The Cable is built in	Reusable Cable is purchased separately from CooperSurgical

(1): If the Patient Return Electrode with the built-in cable is purchased, the Patient Return Electrode Cable is plugged directly into socket 8.

(2): If the Patient Return Electrode and a separate Reusable Cable are purchased, the Patient Return Electrode is connected to the Reusable Cable, which is then plugged into socket 8.

## Section 3 LEEP PRECISION Generator Features

- Electronic power level control
- Digitally derived waveforms
- Low inherent leakage rate
- Patient return monitor disables power automatically if Patient Return Electrode is loose
- Regulated system performance traceable to NIST Standards
- Automatic smoke evacuator signal function sends signal to detector for automatic smoke evacuator actuation when generator is activated
- Self-monitors all detectors

## Classification

<u>Model</u>	<u>Safety Class</u>	<u>Type</u>
LEEP PRECISION Generator	I	BF

- Do not get fluid into the LEEP PRECISION Generator. Should any liquid or solid object fall into the unit, unplug the unit and call Technical Support.
- The LEEP PRECISION Generator is suitable for intermittent operation with a two-minute ON and six-minute OFF duty cycle.
- The LEEP PRECISION Generator is classified as normal equipment (IPX0) according to protection against ingress of water.

# Section 4 Front and Rear Panels of the LEEP PRECISION Generator and Integration Unit

## 4.1 Front Panel of the LEEP PRECISION Generator (Cart not shown)

(Colored, numbered boxes are also located later in this manual.)


Figure 6 (Front Panel)

<p><b>Controls 1</b></p> <p>Operator Pushbuttons</p> <ul style="list-style-type: none"> <li>a) PURGE* Turns LEEP PRECISION Smoke Evacuator on for 5 sec.</li> <li>b) RAISE Increases power setting</li> <li>c) LOWER Decreases power setting</li> <li>d) WAIT Standby mode</li> <li>e) BLEND Cutting mode selection</li> <li>f) CUT Cutting mode selection</li> <li>g) COAG Cutting mode selection</li> </ul> <p>Foot Switch Socket 5</p> <p>A/C Power (rear panel) 6</p> <p>Active Electrode Socket 7</p> <p>Dispersive Electrode Socket 8</p> <p>Remote Pencil Buttons (not shown)</p>	<p><b>Indicator Lights 2</b></p> <p>WAIT = Red CUT = Yellow BLEND = Yellow COAG = Blue</p> <hr/> <p><b>Power Indications 3</b></p> <p>Digital LED Display</p> <hr/> <p><b>Acoustical</b></p> <p>Piezo Annunciator (internal)</p> <hr/> <p><b>Pilot Lights 4</b></p> <p>PATIENT PAD LOOSE R/F LEAKAGE ERROR SYSTEM READY R/F ON</p>
--	--

\* The LEEP PRECISION Smoke Evacuator may also be tested or used without the LEEP PRECISION Generator output by pressing the "Purge Button" on the front panel. The motor will stay on for five (5) seconds after the button is released. Most purge operations (removing excess smoke) take only a few seconds.

## 4.2 Rear Panel of the LEEP PRECISION Generator


Figure 7 (Rear Panel)

### Rear Panel

LEEP PRECISION Generator  
ON/OFF Switch **6**

A/C Power Cord Connector **9**

Fuse Holder **10**

“Communication” Infrared LED between  
the LEEP PRECISION Generator and  
LEEP PRECISION Integration Unit (ON/  
OFF) **11**

### Symbols on the LEEP PRECISION Generator

 Classification I Type BF protected against defibrillator effects

 Floating output circuit (Applied Part)

 Cautions – consult this manual for safety precautions

 Dangerous voltage

 Equipotentiality Symbol (Ground)

## 4.3 LEEP PRECISION Integration Unit


Figure 8

### LEEP PRECISION Integration Unit

Main Power Switch (ON/OFF) **12**

The Power cord connects here **13** and other end goes to the wall outlet

The Connector with Cord **14** is permanently mounted on the LEEP PRECISION Integration Unit; the other end plugs into the LEEP PRECISION Smoke Evacuator

### IMPORTANT

The LEEP PRECISION Integrated System/Generator user must be thoroughly trained in the techniques of Loop Electrosurgical Excision Procedures (LEEP). This system has been designed for use with only the CooperSurgical LEEP PRECISION Electrosurgical Accessories. DO NOT use this equipment for any purpose other than that for which it has been designed. See Warnings and Caution statements throughout this manual.

## Section 5 Professional Use Guide

This manual contains information about the proper procedures for inspecting and preparing the LEEP PRECISION Generator before its use as well as for its care and storage after use.

This manual does not describe how an actual procedure is to be performed, nor is it meant to teach a beginner the proper technique or any of the medical considerations regarding the use of this equipment. CooperSurgical strongly recommends that the prospective user obtain appropriate training before using this equipment, as improper use could be potentially hazardous to the patient and the user.

**This device SHOULD NOT be used without proper training.**

**Training in the use of electrosurgical equipment should include:**

- 1. A review of the published literature regarding the procedure of interest**
- 2. Attendance at a course or courses offered by physicians experienced with the Loop Electrosurgical Excision Procedure**
- 3. Hands-on preceptor training from an experienced practitioner**

**Read this entire manual carefully to become familiar with each of the controls and features before making any attempt to use the equipment clinically.**

Instructions contained in the operating manuals of any equipment to be used in conjunction with this equipment must be followed to avoid any possible hazard from incompatibility.

Failure to thoroughly understand and follow the instructions given in this manual may result in serious injury to the patient and/or the operator. Failure to follow the instructions given in this manual may result in damage to or malfunction of this equipment.

No long-term follow-up studies as to recurrence rates with this device have been performed. The effects of Loop Electrosurgical Excision Procedure on pregnancy outcome are not known.

**SAFETY PRECAUTIONS MUST ALWAYS BE EXERCISED WHEN USING ELECTRICAL EQUIPMENT TO PREVENT OPERATOR/PATIENT SHOCK, FIRE HAZARD, AND EQUIPMENT DAMAGE.**

**CAUTION:** U.S. Federal law restricts this device to sale by or on the order of a physician. This device SHOULD NOT be used without proper training and preceptorship.

If any questions arise regarding the information contained in this manual, the operation, or safety of the equipment or service, please contact CooperSurgical.

### 5.1 LEEP PRECISION Cart Use

**VERY IMPORTANT:** The LEEP PRECISION Integrated System must only be moved by grasping and holding firmly onto the LEEP PRECISION Cart Handles to ensure the Cart does not tip over.

The LEEP PRECISION Cart must only be pushed or pulled using the handles on the front of the LEEP PRECISION Cart to ensure stability during transportation.

### 5.2 Indications

- Cervical conization
- Large Loop excision of the transformation zone (LLETZ) in the diagnosis and treatment of some cervical intraepithelial neoplasias (CIN) and dysplasias
- External anogenital lesion
- Large vaginal intraepithelial neoplastic (VAIN) lesions

### 5.3 Contraindications

- Pregnancy
- Known or suspected cervical changes secondary to DES (diethylstilbestrol) intrauterine exposure
- Acute or active inflammation of the cervix, endometrium, fallopian tube, ovary or peritoneum (cervicitis, endometritis, tubo-ovarian inflammatory disease or pelvic inflammatory disease)
- Invasive cancer that is visible on examination
- Less than three months postpartum

### 5.4 LEEP Procedure and Technique

It is recommended that the patient be provided with a brief description of the procedure and the equipment that will be used.

### 5.5 Safety Precautions

1. The LEEP PRECISION Integrated System should be used only by a thoroughly trained physician in an adequately equipped medical facility.
2. Replacement accessories and Patient Return Electrodes should be kept on hand as defective active accessory or Patient Return Electrodes can result in substandard performance of this equipment.
3. The LEEP PRECISION Integrated System should be connected only to a properly grounded receptacle. NEVER use an adapter that defeats the ground of the built-in three (3) prong plug.
4. Care must be exercised when handling liquids around electrical equipment. DO NOT attempt to operate this equipment if liquids have spilled on any part of the LEEP PRECISION Integrated System. DO NOT use flammable liquids around electrical equipment.
5. This equipment should never be used in conjunction with other equipment for which safety against leakage current has not been established.
6. This equipment should only be used with CooperSurgical LEEP PRECISION Accessories and CooperSurgical LEEP Disposable Accessories.
7. When this equipment is operated:
  - a. A Patient Return Electrode of adequate surface area MUST be properly attached to the patient or the risk of accidental burns will exist.
  - b. The Patient Return Electrode should be placed as close as possible to the site of active accessory use, but must NEVER be placed so as to allow the patient's heart to be in the pathway from the active accessory to the return electrode!
8. The user should thoroughly understand the principles and use of radio frequency (RF) current before using this equipment. This understanding is essential to avoid the hazard of shocks or burns to the user and/or the patient.
9. The instructions for use described in this manual must be followed; otherwise, compromised safety, malfunction, injury to the operator and/or patient, or costly damage to the LEEP PRECISION Integrated System may occur.
10. There are no user-serviceable parts within the housing. Repairs to this equipment should only be performed by authorized CooperSurgical service personnel. For service information, please contact CooperSurgical (refer to Section 19).

## 5.6 Electrosurgical Procedures

This section provides only general information about the use of electrosurgical devices. Only the user can evaluate the clinical factors involved with each patient and determine if the use of this equipment is indicated. The user must then decide on the specific technique and procedure that will accomplish the desired clinical effect.

### WARNING

LEEP PRECISION Generators are designed to allow the controlled destruction of tissue and are inherently dangerous if operated improperly.

REPORTED PROBLEMS DUE TO IMPROPER OPERATION DURING ELECTROSURGICAL PROCEDURES HAVE INCLUDED:

- Inadvertent activation with resultant tissue damage at the wrong site and/or equipment damage
- Alternate current pathways resulting in burns where the patient, physician, or assistant is in contact with exposed metal
- Explosions caused by electrosurgical sparking in a flammable gas mixture (i.e., explosive anesthetic gases and the inappropriate use of alcohol and other flammable liquids)
- Perforation and massive hemorrhage

A proper Patient Return Electrode pathway is extremely important during any monopolar electrosurgical procedure. Every effort must be made to ensure that, throughout the electrosurgical procedure, an adequate surface area is provided and remains in proper contact with the patient to reduce the current density below a level that might cause inadvertent tissue damage where the Patient Return Electrode has been applied.

### 5.6.1 Electrosurgical Tissue Effect

Delivery of continuous sinusoid waveform currents through a small electrode at appropriate power levels can cause rapid heating of the intracellular fluids in the cells in close proximity to the electrode, turning these fluids to steam. The significant increase in volume (approximately five times) causes cellular structures to rupture, creating the clinical effect of CUT, with little or no hemostatic effect along the margin of the divided tissue. Delivery of short-duration pulses of RF currents through a small electrode at appropriate power levels can cause heating of intracellular fluids at a more gradual pace. This allows evaporation of these fluids without rupturing the cellular structure, creating the clinical effect of desiccation, or COAG, without the division of tissue.

By varying the pulse to an intermediate duration, it is possible to get a clinical effect that combines or blends the clinical characteristics of CUT and COAG, yielding the effect referred to as BLEND, in which tissue is divided with a desirable amount of hemostasis along the margins of the divided tissue.

The LEEP PRECISION Generator has output load characteristics that cause the electrosurgical effects to remain consistent throughout the procedure. However, under some circumstances, it may be necessary to readjust the power settings during the procedure.

The electrosurgical effect may vary throughout the procedure, requiring the operator to adjust the LEEP PRECISION Generator relative power settings.

**5.6.2 Select the Output Mode (i.e., CUT, BLEND or COAG by Pushing the Corresponding Button **1** e, f, or g)**

**NOTE:** Numbers in **blue boxes** refer to those component parts in Figure 6.

<b>Output Mode</b>	<b>Waveform Description</b>	<b>General Effect</b>
CUT	Continuous 495 kHz sinusoid with minimal modulation	Cutting without Hemostatis
BLEND	Interrupted 495 kHz sinusoid intermediate-duty cycle	Cutting with minimal Hemostatis
COAG	Bursts of 495 kHz sinusoid short-duty cycle	Coagulation without Cutting

**5.6.3 Set the Output Power level [confirmed on the digital display ( **1** b and c and **3** ) by using the Output Power selector buttons as desired**

**WARNING**

The degree and speed of electrosurgical effect is largely dependent on Current Density at the point of contact of the Active Electrode. Loop Electrosurgical Excision Procedure electrodes from other manufacturers may vary in the diameter, thickness, size, and configuration of the cutting wire. This may result in significant changes in the electrosurgical effect at a given output power level setting. The use of CooperSurgical LEEP PRECISION Electrodes is recommended.

**5.6.4 Attendant Risks of High Frequency Electrosurgery**

Electrosurgery produces sparks that can ignite flammable materials. This includes solvents, adhesives, gauze, cotton, and also liquids and gases.

Further risks include:

- Pooling of flammable agents under the patient in body depressions or in body cavities.
- Oxygen-saturated gauze and cotton wool
- Accumulated oxygen gas: under covers, drapes, or clothing
- Ignition of endogenous gases in or near body cavities or osses

**5.6.5 Electrosurgery**

Electrosurgery is the use of energy created by a high frequency alternating current. The tissue's resistance to the passage of this current creates heat internally in the tissue, as in diathermy. Monopolar electrosurgery involves the use of two electrodes (an Active Electrode and a Dispersive Electrode) of greatly different sizes. This results in greatly increased current density at the point of the smaller electrode. While the electrode itself remains cold, the highly concentrated high frequency energy creates molecular heat inside each cell. By the choice of electrodes and selection and adjustment of the current, the operator controls the effect of this energy on the tissues to achieve the desired results.

**CAUTION**

**FOR MONOPOLAR SYSTEMS, IT IS RECOMMENDED TO KEEP THE VOLTAGE/POWER AS LOW AS POSSIBLE TO ACHIEVE THE DESIRED END EFFECT (due to the potential for capacitive coupling and inadvertent burning at higher voltages).**

With electronically generated electrosurgery current it is necessary to convert the 50/60 Hz alternating current available at the wall outlet to the high frequency current needed for electrosurgery. This conversion is accomplished by a high frequency generator. The high

frequency waveforms are precisely controlled for the various modes of operation. The peak average and RMS values of the desired waveforms are generated in accordance with scientifically and empirically derived standards that are known to result in the desired effects.

### **5.6.6 Fundamentals of Electrosurgery**

As with any technical instrumentation or equipment, there are some fundamental principles that should be learned in order to use electrosurgery safely and effectively. These principles apply generally to all procedures in which electrosurgery is used. Users should familiarize themselves with them.

### **5.6.7 The Active Electrodes**

Be sure the electrodes are securely seated and firmly held in the handpiece so that the entire electrode shank is covered.

In general, straight wire electrodes are used for incisions and for removing fine tissue. Loops are used for removing heavier tissue, planing, and contouring. Ball Electrodes are used for coagulation.

Keep the electrode clean while operating. Tissue shreds and debris on the electrode reduce the effectiveness of the current and, by impeding the passage of the electrode through the tissue, slow down the stroke. This creates unnecessarily heavy coagulation, which can cause sloughing and delay healing.

Inspect the electrodes for proper insulation.

### **THIS DEVICE SHOULD ONLY BE USED BY PHYSICIANS TRAINED IN ELECTROSURGERY.**

The practitioner who lacks experience should not attempt the procedures described in the following pages based solely on this information.

**NOTE:** The best initial effect is accomplished with the cutting wire in only light contact with tissue. Tight pressure may cause desiccation of the tissue and will delay the start of the cutting effect.

### **5.6.8 Thermal Effects on Tissue Treated with Loop Electrodes**

Thermal effects on tissue specimens may include:

- Thermal coagulation injury of the cervix, up to one-third the thickness of normal epithelium of the cervix
- Fragmentation of squamous epithelium of the cervix attributable to long exposure periods along the excision site, which allow heat to dissipate laterally
- Partial coagulation of the endocervical epithelium because of lateral heat radiation. The Loop Electrosurgical Excision Procedure may produce thermal effects at the periphery of the excised tissue and may make histopathologic interpretation difficult or impossible and, therefore, may not allow accurate diagnosis and may obscure the need for further treatment.

## Section 6 Electrosurgical Precautions

The safety and effectiveness of electrosurgery is dependent to a large degree upon the skill of the user/operator. It is important that the user/operator read, understand, and follow the operating instructions supplied with the CooperSurgical LEEP PRECISION Integrated System as well as thoroughly understand the principles and use of electrosurgical systems.

### WARNING

Electrosurgery uses radio-frequency energy to cut and coagulate tissue. Because of the sparking and heat associated with electrosurgery, do not use with flammable anesthetics or other flammable gases, near flammable fluids or objects, or with oxidizing agents.

- Skin-to-skin contact, for instance between the patient's arm and body, should be avoided by the placement of an appropriate separating device, such as two to three inches of dry gauze. This will reduce the potential for alternate site burns.
- If monitoring, stimulation, imaging, or similar devices are used simultaneously with electrosurgery, the monitoring electrodes must be placed as far as possible from the electrosurgery site and the Patient Return Electrode. Position the Patient Return Electrode close to the electrosurgery site, for example, on the thigh when treating the cervix. **NOTE:** Monitoring Needle Electrodes are not recommended.

## WARNINGS

- **THIS DEVICE IS A BF-TYPE DEVICE. BF-TYPE EQUIPMENT IS A B-TYPE DEVICE WITH AN F-TYPE APPLIED PART.**

A B-Type device is a piece of equipment providing a particular degree of protection against electric shock, particularly regarding:

- Allowable leakage current
- Reliability of the protective earth connection.
- An F-Type device is an isolated applied part. The applied part is isolated from all other parts of the equipment to such a degree that the patient leakage current allowable in single fault condition is not exceeded when a voltage equal to 1.1 times the highest rated mains voltage is applied between the applied part and earth.

- **AVOIDING UNINTENTIONAL PATIENT BURNS**

1. Do not activate the electrode until all components are in place.
2. Be certain to have an unobstructed path and view.
3. Use only non-conductive ancillary instruments.
4. Remove patient jewelry, metal drape clips, (etc.) and needles.
5. A Patient Return Electrode must be used. The entire area of the Patient Return Electrode should be reliability attached to the patient's body and as close to the operating field as possible.
6. Place Foot Switch in safe position, out of traffic.
7. Turn console OFF or place in WAIT state when not in use.
8. Do not use unit on small appendages due to current density. Bipolar techniques may be desirable for use on small appendages.

- **THERMAL INJURY AND DEFECTS OF TISSUE TREATED WITH LOOP ELECTRODES**

Possible injury to cervical tissue may include:

1. Thermal coagulation injury of the cervix, up to one-third the thickness of normal cervical epithelium,
2. Fragmentation of squamous cervical epithelium attributable to long exposure periods along the excision site that allows heat to dissipate laterally, and
3. Partial coagulation of the endocervical epithelium because of lateral heat radiation.

Loop Electrosurgical Excision Procedures may also produce thermal defects at the periphery of the excised tissue and may make histopathologic interpretation difficult or impossible, preventing accurate diagnosis and the evaluation of the need for further treatment.

- **FLAMMABLE GASES / SOLUTIONS**

1. When using a flammable preparatory solution, such as alcohol, be certain that the solution has completely evaporated before using electrosurgery.
2. Electrosurgery should not be used in the presence of flammable or explosive gases.
3. It is recommended that only non-flammable agents be used for cleaning and disinfection.

SMOKE FROM PROCEDURES CONTAINS LIVE CELLS THAT MAY PRESENT BIOHAZARDS (HIV, ETC.); USE A SMOKE EVACUATION DEVICE.

CONSULT THE PACEMAKER MANUFACTURER'S LITERATURE AND THE PATIENT'S PHYSICIAN BEFORE USING THIS INSTRUMENT ON A PATIENT USING A PACEMAKER.

### CAUTIONS

- U.S. Federal law restricts this device to sale by or on the order of a physician. This device SHOULD NOT be used without proper training and preceptorship.
- For monopolar systems it is recommended to keep the voltage/power as low as possible to achieve the desired end effect (due to the potential for capacitive coupling and inadvertent burning at high voltages).
- Adequate anesthesia is indicated for all electrosurgical procedures.
- Inspect loop and ball for adequate insulation.
- Avoid using an electric extension cord (danger of separation).
- Use a Smoke Evacuator with this device.
- This device generates high-frequency energy that could, under some circumstances, cause interference to other equipment.
- For situations where the MAXIMUM OUTPUT VOLTAGE is less than or equal to 1,600 V, ASSOCIATED EQUIPMENT and ACTIVE ACCESSORIES should be selected that have a RATED ACCESSORY VOLTAGE equal to or greater than the MAXIMUM OUTPUT VOLTAGE.
- ASSOCIATED EQUIPMENT and ACTIVE accessories should be selected with RATED ACCESSORY VOLTAGE  $\geq$  MAXIMUM OUTPUT VOLTAGE when the smaller of variable  $y$  [see below] or the number 6 is  $\leq$  CREST FACTOR for that HF SURGICAL MODE.
- When MAXIMUM OUTPUT VOLTAGE ( $U_{max}$ ) is  $> 1,600$  V, and the CREST FACTOR is  $<$  the variable  $y$  calculated below, indicating that any ASSOCIATED EQUIPMENT and ACTIVE ACCESSORIES used in such mode or setting must be rated to withstand the combination of actual voltage and CREST FACTOR.

$$y = \frac{U_{max} - 400 \text{ [Volts]}}{600 \text{ [Volts]}}$$

- No modification of this equipment is allowed.
- No customer-serviceable parts.
- Only replacement of power cords or fuses can be done by the user. Only replace with these items with the exact replacement part(s) available from CooperSurgical.

## Section 7 LEEP PRECISION Patient Return Electrode

A Patient Return Electrode must be used to operate the LEEP PRECISION Integrated System/Generator. The Patient Return Electrode must be of the Split or Dual type for the safety features to prevent R/F burns due to poor adhesion.

### NOTE

For optimum patient safety with the use of the LEEP PRECISION Integrated System/Generator, it is strongly recommended to use ONLY authorized and genuine CooperSurgical LEEP PRECISION Split/Dual Patient Return Electrodes. No other Patient Return Electrode has been tested for use with this system.

### WARNING

**Patient injury may result. Do not use a single/not split Patient Return Electrode as it is not continuously monitored for integrity.**

## Section 8 Turning on the LEEP PRECISION Integrated System

**NOTE:** Numbers in **blue boxes** refer to those component parts in Figures 6 and 8.

### 8.1 Practice

The following section describes cutting and coagulating. Practice methods are first suggested, including how to adjust the settings for optimal cutting and coagulating. Then general techniques are described, which should help determine and develop the best specific techniques.

1. Refer to Section 2 for the general setup of all parts of the LEEP PRECISION Integrated System and make sure every part of the system is connected together as described.
2. Connect the A/C power cord(s) accordingly. If not already plugged in, plug the removable Foot Switch into the Front Panel (at **5**), if foot control is desired. Set it up in a convenient, low-traffic area. Refer to Figures 6 and 8.
3. Make sure the ON/OFF switch **12** on the LEEP PRECISION Integration Unit is in the OFF position. Plug the LEEP PRECISION Integrated System into a grounded 120 VAC or a 230 VAC electrical receptacle. (The Power Cord receptacle **13** is shown in Figure 8.) This is an important safety feature.
4. Use only approved CooperSurgical LEEP PRECISION disposable accessories.
5. Remove the disposable Patient Return Electrode from its packaging and attach it to the patient in accordance with the instructions. Be sure to follow all instructions regarding the location of the Patient Return Electrode to ensure that good contact is made with the patient.
6. Attach the Patient Return Electrode securely to the Dispersive Electrode Socket **8** on the LEEP PRECISION Generator's front panel.
7. Insert the Active Electrode into the Active Electrode Socket **7** on the top right side of the LEEP PRECISION Generator's front panel. Either the Hand Switch-Operated Electrode (Pencil) or the Foot Switch-Operated Electrode (Pencil) may be used.
8. Insert the selected Electrode (Needle, Loop, Square or Ball) into the Active Electrode.
9. Set the ON/OFF switch **12** to ON. A red indicator lamp on the console (R/F ON) **4** will indicate that the hold mode has been initialized. After Self-Check, the SYSTEM READY LED **4** (green) will illuminate.

#### NOTES

If the LEEP PRECISION Integrated System does not cycle to the SYSTEM READY status and a series of prolonged audible beeps are heard, turn the system off and check for faulty hand, foot, or panel switches. See the **TROUBLESHOOTING** section (Section 12) of this manual for additional information.

**Do not be alarmed if the LEEP PRECISION Smoke Evacuator turns on at this point. This is the unit cycling.**

10. Ensure patient is connected securely to the Patient Return Electrode before beginning.

#### WARNING

**POSITIVE CONTACT MUST BE COMPLETE BETWEEN THE PATIENT AND THE PATIENT RETURN ELECTRODE. IF NOT, SEVERE BURNING OF TISSUE MAY RESULT. THE ENTIRE AREA OF THE PATIENT RETURN ELECTRODE SHOULD BE RELIABLY ATTACHED TO THE PATIENT'S BODY AND AS CLOSE TO THE OPERATING FIELD AS POSSIBLE.**

11. The output is selected by the mode selector buttons. When a mode is selected, the intensity should be set by depressing the RAISE (**1**b) or LOWER (**1**c) button. If either CUT (pure cut/blended cut) or COAG mode is selected and the Foot Switch is depressed, the output will be active. When the Foot Switch or Active Electrode is actuated, the R/F ON lamp **4** will illuminate and the beeper will emit a pulsing tone.

12. The LEEP PRECISION Integrated System is ready for use and the SYSTEM READY LED **4** (green) will illuminate.

If the operator feels the power output is not sufficient, the Patient Return Electrode should be checked to see that good contact is made with the patient before increasing the power setting.

13. **CROSS CONTROL:** The sequence of operation of the LEEP PRECISION Integrated System is nonpreferential; that is, either the Foot Switch (pedal) or the Hand Switch will activate the output (providing all safety measures are met). The Foot Switch and the Hand Switch cannot be operated simultaneously.

#### IMPORTANT

If using the Foot Switch, the button on the front panel must be used to control the OUTPUT mode. When the COAG mode (**1**g) is selected on the front panel, the CUT or BLEND modes cannot be invoked by the Hand Switch (yellow) Button. The CUT mode, if desired, must be selected on the front panel and then invoked either by the Hand Switch- or the Foot Switch-Operated Pencil.

#### WARNING

Positive contact must be complete between the patient and the Patient Return Electrode. If a Split/Dual Patient Return Electrode is not securely affixed to the patient, the PATIENT PAD LOOSE indicator LED **4** (red) will indicate an improper return path. The PATIENT PAD LOOSE safety circuitry disables the output when an unsafe condition occurs in the Patient Return Electrode circuit. Poor connections or insufficient contact area can cause R/F burns. The entire area of the Patient Return Electrode should be reliably attached to the patient's body and as close to the operating field as possible. Use only approved CooperSurgical LEEP PRECISION Dual/Split Patient Return Electrodes. Do not use a single (not split) Patient Return Electrode as it may cause severe tissue burns in the vicinity of the Patient Return Electrode.

14. **CABLE PLACEMENT:** High frequency cables should be routed clear of other instruments and arranged for minimum length. Cables to surgical electrodes should be positioned to prevent contact with the patient or other leads. Temporarily unused active electrodes should be stored in a location isolated from the patient.

#### NOTE

The RF nature of the electro surgical generator has the potential to electrically couple to the ground and therefore could introduce RF current leakage (R/F LEAKAGE section). Attention should be given to the routing of the wires of the system accessories to prevent this situation. Here are some of the suggestions that need to be considered to address this issue:

- The wires should be hanged (i.e. Cart, IV pole, Surgical Bed, etc.) to something that is not grounded (coated and not bare metal).
- The wires should be straight and not coiled.
- The wires can be routed together to increase the effectiveness of the cable routing.

15. **R/F LEAKAGE:** The LEEP PRECISION Integrated System/Generator has a sensitive error detection mode that disables the output and indicates that an error condition exists by illuminating the front panel R/F LEAKAGE ERROR LED **4** (red). This LED will indicate that an undesirable path exists from the Patient Return Electrode or Active Electrode to ground. In this context, ground means any undesirable path or place. If such a path exists, the patient should be inspected for any contact with metallic paths to tables, chairs, or other non-isolated objects. The purpose of this feature is to prevent R/F current from flowing in undesirable places where it could either cause burns or damage to ancillary equipment.
16. **PATIENT CONTACT:** Direct operator-to-patient skin-to-skin contact could create an undesired leakage path. Avoid skin-to-skin contact by padding with dry gauze as necessary. The operator should wear electrically insulating gloves. Unless absolutely necessary, the patient should not be touched when R/F power is activated. Patient should not come into contact with grounded metal parts or parts with a higher capacitance to ground (i.e., operating table supports). Use of anti-static sheeting is recommended for this purpose.

## 8.2 The Power Setting

Select the operating mode (CUT, BLEND, COAG, or WAIT) by using the mode selector buttons. When a mode is selected, the intensity should be set by depressing the POWER RAISE ( **1** b) or LOWER ( **1** c) button that controls the digital display. When a mode is selected, the corresponding LED illuminates in the corner of each button to verify which mode is engaged.

To ensure correct operation of the instrument, the colors are correlated to international standards:

**WAIT = Red**  
**BLEND = Yellow**  
**CUT = Yellow**  
**COAG = Blue**

When the Foot Switch is depressed, or the Hand Switch Buttons are pressed, the R/F ON **4** (yellow) will illuminate. If either the CUT ( **1** f) (pure or blend) or COAG ( **1** g) mode is selected and the Foot Switch is depressed, the output will be active. In the WAIT ( **1** d) mode there will be no output. The Piezo Electric Annunciator will indicate with an intermittent tone when the R/F Power is present at the output. Once contact is made with the electrode to the tissue and a small current is flowing, the tone will become steady.

The ideal power setting is the minimum setting at which cutting and coagulating meet the criteria described in Sections 8.4 – 8.6. If the power is too high the tissue will be discolored and there will be considerable sparking when the electrode makes contact with the surface. If the power is too low the electrode will drag through the tissue, tearing and burning instead of cutting cleanly, and will pick up shredded and torn tissue.

Power requirements will vary with the type and size of the electrode, the area of electrode surface in contact with the tissue, the nature of the tissue, whether cutting or coagulating, and the depth of the incision desired. Larger electrodes, deeper incisions, and tough fibrotic tissue are some indications for higher power settings.

### CAUTIONS

Do not use fine Needle Electrodes at high settings as they may be severely damaged by such use. When practicing with these electrodes, work up gradually from a low setting. Generally, you should not exceed a setting of 40 with fine Needle Electrodes.

Do not use monitoring Needle Electrodes during electrosurgery. If monitoring electrode cables are used, they should be placed as far away from electrosurgical cables AS POSSIBLE. If it is absolutely necessary to cross over any monitoring cable or patient-connected tubes, such as IV lines, the crossing should be at right angles. Use of monitoring systems incorporating high frequency current limiting devices is recommended.

## 8.3 Cutting Techniques

When cutting, ALWAYS activate the electrode by stepping on the Foot Switch or pressing the yellow button on the Active Electrode (Pencil) BEFORE making contact with the tissue.

Plan the stroke. Before activating the electrode, take one or two practice strokes to be sure you can complete the planned stroke comfortably and correctly. At this time you can evaluate the size and shape of the electrode and the speed and depth of the stroke.

When the practice stroke is comfortable, step on the Foot Switch or press the yellow button on the Active Electrode (Pencil) and make the planned cut.

Use a smooth brushing motion without pressure. The electrode should pass through the tissue, WITHOUT DRAGGING, at a deliberate but not slow speed.

Keep the electrode moving. Prolonged contact with any one part of the tissue can create excessive coagulation.

### 8.3.1 Cutting

**Blend:** Cutting is done not by the electrode, but by the high frequency energy concentrated at the electrode's tip. This high frequency energy generates molecular heat in each cell to the point at which the fluids in the cell vaporize and the cell explodes. By applying this energy to individual cells in sequence, that is, by moving the electrode continuously through the tissue, the line of destruction is limited and the cutting effect is realized. At the same time, the capillaries are sealed, resulting in almost bloodless cutting, hence the term "blended cut."

**Pure cut:** This is almost like using a cold scalpel with very little or no hemostasis.

## 8.4 Criteria of a Good Cutting Technique

There are three criteria of a good cutting technique:

1. The electrode should "float" through the tissue without dragging or resistance.
2. There should be only very slight, if any, change in the tissue color due to dehydration or charring.
3. No tissue shreds should adhere to the electrode.

## 8.5 Coagulating

The LEEP PRECISION Generator is designed for high fulguration so tissue bleeders can be sealed without burning the uninvolved tissue. Coagulation takes place when the high frequency current is applied to the tissue with a current density sufficiently concentrated to dehydrate the cells and coagulate their organic contents, but without penetrating deeply into the tissue.

This procedure is almost self-limiting because the surface coagulation first created protects the underlying tissue against excessive depth of coagulation.

Coagulation appears as a white spot on the surface of the tissue, emanating from the point of contact to the Ball Electrode. The depth of coagulation is approximately equal to the lateral spread of coagulation.

In contrast to cutting, when coagulating, the Ball Electrode should make contact with the tissue before the physician depresses the Foot Switch or Hand Switch.

Start with a LOW power setting (15), changing the mode button to COAG.

Make light contact with the surface of the tissue with the Ball Electrode. Depress the Foot Switch or Hand Switch to activate the current, coagulating the area for several seconds, then release the Foot Switch or Hand Switch and remove the electrode from the tissue.

Note the appearance of the tissue, which should appear blanched. Turn the power intensity up step by step, repeating the procedure. Observe the characteristics of the coagulated tissue with each setting. The degree of coagulation obtained at a particular setting will vary with different tissue and different conditions. Determine which setting is best for various procedures. When treating any condition, if proper coagulation is not apparent after one application, immediately increase the intensity setting as long as the Patient Return Electrode is properly in place.

### CAUTION

Never repeatedly apply the current to the same area; this may create heat and seriously damage underlying tissues.

## 8.6 Coagulating Technique

The Ball Electrode can control hemorrhage. When applied for a second or two, electronic coagulation current will coagulate small capillaries easily. Larger vessels can be picked up with hemostats, which in turn may be touched by the Ball Electrode using electronically generated current. This will seal all vessels that are held in the tips of the forceps.

### NOTE

Always make contact with the tissue before depressing the Foot Switch or Hand Switch when using the coagulation current.

## 8.7 Technique Guidance

### 8.7.1 For OB/GYN Procedures

1. The endocervix may not be included in the loop excision, and the results of endocervical curettage (ECC) may not be predictive of either residual or invasive disease after loop excision procedures. If the ECC is positive for dysplasia, a standard cone biopsy should be considered.
2. Loop Electrosurgical Excision Procedures performed with small diameter wire Loop Electrodes produce multiple small pieces of cervical tissue and may provide a less acceptable tissue specimen for histopathologic analysis.
3. Larger lesions involving multiple quadrants of the cervix are more difficult to remove with Loop Electrodes.

### 8.7.2 For Severed Vessels

1. Clamp the bleeder with a hemostat.
2. Using electronically generated current, touch a Ball Electrode to any part of the hemostat.
3. Depress the Foot Switch Pedal to activate the current for several seconds, then release the Pedal and remove the electrode. Be sure the current is set as high as is necessary to achieve coagulation.
4. After the application of coagulating current, remove the hemostat. Bleeding should have stopped; if not, repeat the procedure. Suture ligation may be required if R/F is unsuccessful.
5. When using electronically generated current, it is safe to hold the hemostat in a gloved hand while holding the electrode handle in the other during this procedure.

### 8.7.3 For Anesthesia

Adequate anesthesia is indicated for all electrosurgical procedures. It is usually advisable to anesthetize tissue adjacent to the intended operative site in the event it becomes necessary to extend the operative area.

### 8.7.4 For Biopsy

The use of electrosurgery for cervical biopsy has advantages in that it seals the capillaries and lymphatics as it cuts. Whenever possible, the specimen should include two to three millimeters of attached normal tissue.

Small (up to one half inch) masses should be removed in one piece. Using a cutting current at a relatively high setting and a Needle Electrode, incise all around the mass in an elliptical pattern, including two to three millimeters of attached normal tissue.

Small masses may also be removed by using a suitable Loop Electrode – large enough so that it can excise the mass as well as two to three millimeters of attached normal tissue at the same time.

Specimens from larger masses should be taken in the form of wedges. Using a Needle Electrode and cutting current, start at the apex of the wedge at the center of the suspicious mass. The base of the wedge should include two to three millimeters of attached normal tissue.

## Section 9 Powering Down the LEEP PRECISION Integrated System

1. At the completion of each procedure, activate the system to ensure safe particle containment. Using gloves and a mask, remove the Pre-Filter, the Reducer, and the used section of Suction Tubing and discard into an infectious waste receptacle (see Cautions). The CooperSurgical Smoke Evacuator should be stored with a new Pre-Filter and Reducer in place on the ULPA filter.
2. Turn off the Power Switch (on the LEEP PRECISION Integration Unit) and unplug from the wall outlet.
3. Dispose of the used sterile Active Electrode (Pencil) and Patient Return Electrode (Dispersive Electrode).
4. Do not throw out the Adapter that is used for the Foot Switch-Operated Electrode. This Adapter is needed each time a new Foot Switch-Operated Electrode is used for each subsequent procedure.
5. If the Reusable Cable was used for the Patient Return Electrode, retain it for each subsequent use of a new Patient Return Electrode.

## Section 10 Maintenance

The LEEP PRECISION Integrated System must be cleaned and disinfected after each use. See Section 14.6 for cleaning and disinfection instructions.

## Section 11 Accessories

Use only CooperSurgical LEEP PRECISION accessories for optimal system performance and patient safety. The list includes, but is not limited to:

- LEEP Loop Electrodes
- LEEP Square Electrodes
- LEEP Needle Electrodes
- LEEP Ball Electrodes
- Patient Return Electrodes
- Hand Switch-Operated Electrodes (Pencils)
- Foot Switch-Operated Electrodes (Pencils)
- 4 mm Reusable Electrosurgical Adapter (for Hand Switch-Operated Electrodes)
- Foot Pedal

The use of non-CooperSurgical-authorized accessories is not recommended as they have not been tested for use with this system.

### 11.1 Disposable Hand Switch Pencils

All CooperSurgical Disposable Hand Switch Pencils are sold sterile and are for single patient use only.

**For the LEEP PRECISION Generator, use only CooperSurgical genuine Disposable Hand Switch Pencils (activated accessory with Finger Switch), which comply with IEC 60601-2-2: 4th edition.**

**Genuine CooperSurgical Disposable Foot Switch-Operated Pencils and electrodes are available, sold sterile, and are for single patient use only.**

If you have any questions or require specific accessories, contact your CooperSurgical Representative. A variety of CooperSurgical disposable electrodes and accessories are available. Use only cables approved by CooperSurgical. **Regular inspection is required for accessories, including electrode cables, for damage to the insulation. If damage is found, the accessory should be replaced to assure safe operation.**

**Only use accessories specifically made for the CooperSurgical LEEP PRECISION Integrated System/Generator that are specifically designed for the safe and proper operation of this unit. Use of alternate accessories or parts is not recommended as they have not been tested and verified and use may result in unsafe operation of this unit.**

## Section 12 Troubleshooting

Problem	Probable Cause	Corrective Action
<p style="text-align: center;"><b>(A)</b></p> <p>ON/OFF switch, when in position, does not illuminate</p>	<ol style="list-style-type: none"> <li>1. Unit not plugged in</li> <li>2. Blown fuse</li> </ol>	<ol style="list-style-type: none"> <li>1. Plug into wall outlet</li> <li>2. Replace fuse</li> </ol>
<p style="text-align: center;"><b>(B)</b></p> <p>PATIENT PAD LOOSE, intermittent power output, or low power</p>	<ol style="list-style-type: none"> <li>1. Faulty application or poor contact</li> <li>2. Improper placement</li> <li>3. Worn insulation on electrical leads</li> </ol>	<ol style="list-style-type: none"> <li>1. Check full contact to patient</li> <li>2. Place in accordance with instructions</li> <li>3. Replace</li> </ol>
<p style="text-align: center;"><b>(C)</b></p> <p>R/F LEAKAGE ERROR</p>	<ol style="list-style-type: none"> <li>1. Electrode touching grounding metal</li> <li>2. Directly grounded</li> <li>3. Insulation missing</li> </ol>	<ol style="list-style-type: none"> <li>1. Remove from metal</li> <li>2. Locate faulty insulation between table/equipment to patient</li> <li>3. Locate faulty insulation</li> </ol>
<p style="text-align: center;"><b>(D)</b></p> <p>SYSTEM READY light (green) not on</p>	<ol style="list-style-type: none"> <li>1. Improper connection</li> <li>2. One or more safety problems as indicated by other lamps or annunciator</li> <li>3. Product failure</li> </ol>	<ol style="list-style-type: none"> <li>1. Connect correctly</li> <li>2. Follow safety instructions</li> <li>3. Return to CooperSurgical</li> </ol>
<p style="text-align: center;"><b>(E)</b></p> <p>R/F ON doesn't light up:</p> <ul style="list-style-type: none"> <li>• No power when Foot Switch is depressed</li> <li>• No power when the switch on the Active Electrode (Pencil) is depressed</li> </ul>	<ol style="list-style-type: none"> <li>1. Electrical current not flowing</li> <li>2. Defective Foot Switch</li> <li>3. Defective Pencil</li> </ol>	<ol style="list-style-type: none"> <li>1. Check all connections</li> <li>2. Replace</li> <li>3. Replace</li> </ol>
<p style="text-align: center;"><b>(F)</b></p> <p>Digital readout does not illuminate</p>	<ol style="list-style-type: none"> <li>1. Unit not connected</li> <li>2. Digital readout circuitry failure</li> </ol>	<ol style="list-style-type: none"> <li>1. Plug into wall outlet</li> <li>2. Return to CooperSurgical</li> </ol>
<p style="text-align: center;"><b>(G)</b></p> <p>System beeps long intervals</p>	<ol style="list-style-type: none"> <li>1. Stuck Panel Switch</li> <li>2. Stuck Foot Switch</li> <li>3. Stuck Hand Switch Button</li> <li>4. A Foot Switch or Hand Switch is active upon power up</li> </ol>	<ol style="list-style-type: none"> <li>1. Turn unit off, press each switch, and make sure they depress freely</li> </ol>

## Section 13 Specifications

### 13.1 LEEP PRECISION Integrated System

#### PHYSICAL

**Size (W x D x H):** . . . . . 22.44 inches x 22.07 inches x 41.22 inches (57.0 cm x 56.06 cm x 104.70 cm)

**Weight:** . . . . . approx. 88 pounds (40 kg) (Cart with integration unit + Generator + Smoke Evacuator)

#### ELECTRICAL

<b>Input Voltage</b>	120 VAC	230 VAC
<b>VA Ratings</b>	1300 VA	1300 VA
<b>Fuses</b>	10A, 250V, T-Type, Slow Blow, (Littelfuse® 0218010 or equivalent)	6.3A, 250V, T-Type, Slow Blow, (Littelfuse® 021806.3 or equivalent)

#### HIGH FREQUENCY OUTPUT\*:

Rated output power is 100 watts, 495 kHz ± 5% into 500 ohms, monopolar only.

Mode	Minimum & Maximum Power Settings	Crest Factor	Maximum Open Circuit Voltage
CUT	10 to 100 watts	1.2	600
BLEND	10 to 100 watts	2.5	600
COAG	10 to 100 watts	6.5	4000

\*Accuracy of output control setting:

For output powers in excess of 10 Watts, the actual power as a function of the load resistance and output control setting shall not deviate from that shown in the diagrams shown later by more than ±20 percent.

**RF Leakage:** 100 MA MAXIMUM AT MAXIMUM POWER SETTINGS

### 13.2 LEEP PRECISION Generator

#### PHYSICAL

**Size (W x D x H):** . . . . . 9.38 inches x 11.63 inches x 4.25 inches (23.8 cm x 29.5 cm x 10.8 cm)

**Weight:** . . . . . 10 pounds, 9 ounces (5.25 kg)

#### ELECTRICAL

<b>Input Voltage</b>	120 VAC	230 VAC
<b>Maximum Current</b>	3.15 amps	1.55 amps
<b>VA Rating</b>	132 VA	242 VA
<b>Power Line leakage</b>	Less than 50 micro-amps	Less than 50 micro-amps
<b>Fuses (5 x 20 mm)</b>	Two T 3.15A 250V Slow-Blow (Littelfuse® 02183.15 or equivalent)	Two T 2.5A 250V Super Slow-Blow (Bussmann® GMD - 2.5A or equivalent)

#### HIGH FREQUENCY OUTPUT\*:

Rated output power is 100 watts, 495 kHz ± 5% into 500 ohms, monopolar only.

Mode	Minimum & Maximum Power Settings	Crest Factor	Maximum Open Circuit Voltage
CUT	10 to 100 watts	1.2	600
BLEND	10 to 100 watts	2.5	600
COAG	10 to 100 watts	6.5	4000

\*Accuracy of output control setting:

For output powers in excess of 10 watts, the actual power as a function of the load resistance and output control setting shall not deviate from that shown in the diagrams shown after the specifications section by more than ±20 percent.

**RF Leakage:** 100 MA MAXIMUM AT MAXIMUM POWER SETTINGS

### 13.3 LEEP PRECISION Smoke Evacuator

<b>Input Voltage</b>	120 VAC, 60 Hz	230 VAC, 50/60 Hz
<b>Fuses</b>	15 amps, 250V, T-Type	10 amps, 250V, T-Type

**DIMENSIONS (W x D x H):** . . . . . 9 inches x 9 inches x 22 inches  
(22.86 cm x 22.86 cm x 53.24 cm)

**WEIGHT:** . . . . . 21 pounds (9.53 kg)

**AIR FLOW:** . . . . . 35 CFM minimum at maximum setting

**FILTER:** . . . . . Pre-Filter, charcoal filter for odor removal and hydrophobic  
ULPA filter rated at 99.999 percent efficiency level for  
0.014 micron particles based on Standard CNC Evaluations

### 13.4 LEEP PRECISION Cart with Integration Unit

**DIMENSIONS (W x D x H):** . . . . . 22.44 inches x 22.07 inches x 41.22 inches  
(57.0 cm x 56.06 cm x 104.70 cm)

**WEIGHT:** . . . . . approx. 54 pounds (24.50 kg)

**See the Power Output Characteristics and Power Output at Various Load Resistance graphs on the next two pages**

### LEEP PRECISION Integrated System Cautions

- No modification of this equipment is allowed.
- No customer serviceable parts
- Only replacement of power cords or fuses can be done by the user. Only replace with these items with the exact replacement part(s) available from CooperSurgical.
- Position the LEEP PRECISION Integrated System to maintain access to the power outlet it is plugged into to assure ease of disconnection from the power source if needed.
- To remove all power from the LEEP PRECISION Integrated System, disconnect the power cord from the wall outlet.
- To avoid risk of electrical shock, this equipment must only be connected to a supply mains with protected earth.

### Environment Conditions




	<u>Use</u>	<u>Shipping and Storage</u>
<b>Environmental Temperature</b>	between +10 °C and +40 °C	between +10 °C and +40 °C
<b>Relative Humidity</b>	between 10% and 90%	between 10% and 90%
<b>Air Pressure</b>	between 700 hPa and 1060 hPa	between 700 hPa and 1060 hPa

### WEEE Directive




In order to preserve, protect and improve the quality of the environment, protect human health and utilize natural resources prudently and rationally – do not dispose of waste electrical or electronic equipment (WEEE) as unsorted municipal waste. Contact local WEEE disposal sites.

### 13.5 Power Output Characteristics

### 13.6 Power Output at Various Load Resistance

## Section 14 LEEP PRECISION Smoke Evacuator

### 14.1 System Description

The CooperSurgical LEEP PRECISION Smoke Evacuator (three-stage air filtration system) is used to remove airborne particulate plumes produced during office and surgical procedures and has the following features:

- Low noise level.
- Triple filtration of air provides efficiency level for 0.014 microns rated at 99.999 percent. This includes a Pre-Filter, a charcoal filter for odor removal, and a final safety filter placed after the charcoal filter.
- High air flow for plume collection.
- Attaches to the CooperSurgical LEEP PRECISION Integrated System.

### 14.2 Initial System Setup

The system arrives in one carton that contains the Vacuum Unit. The disposable accessories necessary to conduct procedures may be purchased separately from CooperSurgical.

### 14.3 ULPA Filter Installation

Tilt the LEEP PRECISION Smoke Evacuator forward and insert the large ULPA Filter cylinder with the Air Flow arrow pointing downward (see Photo A).

Photo A

### 14.4 Operating Instructions

Photo B

#### 14.4.1 Installing the Pre-Filter

Prepare the unit for use by inserting a clean, disposable Pre-Filter (Single Patient Use) onto the ULPA Filter cylinder (see Photo B). Be sure this device is firmly seated.

#### 14.4.2 Connecting the Tubing

There are two tubing hookup options:

1. For procedures requiring close-proximity plume removal (i.e., vaginal speculum)
2. For procedures requiring open-area plume removal (i.e., external lesions)

#### For procedures requiring close proximity plume removal (i.e., vaginal speculum)

Assemble the  $\frac{3}{8}$ -inch Reducer (Catalog number 6083) onto the port on the disposable Pre-Filter top with a slight twisting motion. Attach one end of an appropriate length of  $\frac{3}{8}$ -inch inner diameter (ID) Evacuation Tubing (Catalog number 6084) to the Reducer Connector and direct the other end to the patient and any appropriate device being used, such as a vaginal speculum equipped with a Smoke Evacuation Adapter (see Photo C).

Photo C

#### For procedures requiring open area plume removal (i.e., external lesions)

1. Assemble the sterile Disposable Evacuation Tubing (Catalog number 6084) directly into the top of the Pre-Filter. Position the opposite end over the site to be treated.
2. Make sure the LEEP PRECISION Smoke Evacuator is plugged into the LEEP PRECISION Integration Unit (refer to Section 2).
3. Push the POWER switch on the LEEP PRECISION Integration Unit to the ON position to start the LEEP PRECISION Smoke Evacuator. The LEEP PRECISION Smoke Evacuator will activate upon use of the Foot Switch or Active Electrode (Pencil).

4. At the completion of each procedure activate the system to ensure safe particle containment. Using gloves and a mask, remove the Pre-Filter, the Reducer, and the used section of Suction Tubing and discard into an infectious waste receptacle (see Cautions). The CooperSurgical Smoke Evacuator should be stored with a new Pre-Filter and Reducer in place on the ULPA filter.

**NOTE:** The expected life of the ULPA filter is 3–6 months depending on usage. It must be discarded into an infectious waste receptacle if plume odor is detected or the suction is diminished.

## 14.5 Maintenance Procedures

No maintenance is required during normal operation other than making sure that ample space is maintained around the unit to allow for free air flow and adequate cooling.

## 14.6 Cleaning and Disinfection

After each patient use, follow the combined cleaning/disinfection instructions provided below for the LEEP PRECISION Smoke Evacuator, the LEEP PRECISION Generator, and the LEEP PRECISION Cart.

### CLEANING

1. Prepare an enzymatic, neutral pH cleaner solution according to the Manufacturer's instructions.
2. Immerse a clean, lint-free wipe into cleaning solution and thoroughly wring. Thoroughly wipe the front panel surfaces in a circular motion. Discard the wipe.
3. Immerse a fresh, clean, lint-free wipe under warm flowing, utility water and thoroughly wring excess water. Thoroughly wipe the front panel surface and crevices for at least 30 seconds. Discard the wipe.
4. Rinse a fresh, clean, lint-free wipe under warm flowing, utility water and thoroughly wring excess water. Thoroughly wipe the front panel surface and crevices for at least 30 seconds. Discard the wipe.
5. Inspect the device for visible soil or debris. If visible soil remains, repeat cleaning steps 2-4 until device is visually clean.

### DISINFECTION

6. Use a 70% IPA wipe, or a low-lint wipe, saturated with 70% IPA, to thoroughly wipe the device. Pay particular attention to difficult to reach areas such as device seams and area between buttons.
7. With a fresh wipe saturated with 70% IPA, thoroughly wipe the device. Allow the wiped and wetted device to stand for a minimum of the disinfectant product manufacture's labeled contact time. Re-wet the wipes as necessary to ensure that all surfaces remain wet for the entire time.
8. Allow the device to air dry. Once dry, visually inspect the device for any residual soil.
9. Place the disinfected device in an appropriate dry storage area.

## 14.7 Cautions

- U.S. Federal law restricts this device to sale by or on the order of a physician.
- This device produces a strong vacuum force and care should be exercised to ensure that the suction control and the position of the inlet end of the suction tubing are adjusted properly to prevent injury to the patient or inadvertent damage to surgical materials.
- The materials in the plume removed by this device are potentially hazardous. Handle according to 29 CFR 1910.1030 and OSHA 3127.1992 (Occupational Exposure to Bloodborne Pathogens) guidelines.
- To prevent a fire or explosion hazard, do not use the system in the presence of flammable or potentially flammable materials.
- Do not allow fluid to be pulled into the system.
- To prevent premature failure of the ULPA filter cylinder, do not operate this device without a disposable Pre-Filter in place.

## Section 15 EMC Compliance Information for the LEEP PRECISION Integrated System

- Medical electrical equipment needs special precautions regarding EMC and must to be installed and put into service according to the EMC information provided in the accompanying documents.
- Portable and mobile RF communications equipment can affect medical electrical equipment.

### GUIDANCE AND MANUFACTURER’S DECLARATION – ELECTROMAGNETIC EMISSIONS

The CooperSurgical LEEP PRECISION Integrated System is intended for use in the electromagnetic environment specified below. The customer or the end user of the CooperSurgical LEEP PRECISION Integrated System should assure that it is used in such an environment.

Emissions Test	Compliance	Electromagnetic Environment Guidance
<p align="center"><b>RF Emissions</b> <b>CISPR 11</b></p>	<p align="center">Group 1</p>	<p>CooperSurgical LEEP PRECISION Integrated Systems use RF energy only for their internal function. Therefore, RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.</p>
<p align="center"><b>RF Emissions</b> <b>CISPR 11</b></p>	<p align="center">Class A</p>	<p>CooperSurgical LEEP PRECISION Integrated Systems are suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.</p>
<p align="center"><b>Harmonic Emissions</b> <b>IEC 61000-3-2</b></p>	<p align="center">Class A</p>	
<p align="center"><b>Voltage Fluctuations / Flicker Emissions</b> <b>IEC 61000-3-3</b></p>	<p align="center">Complies</p>	

## Section 16 Guidance and Manufacturer's Declaration – Electromagnetic Immunity

The CooperSurgical LEEP PRECISION Integrated Systems are intended for use in the electromagnetic environment specified below. The customer or the end user of the CooperSurgical LEEP PRECISION Integrated System should assure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environmental Guidance
Electromagnetic discharge (ESD)  IEC 61000-4-2	±6 kV contact  ±8 kV air	±6 kV contact  ±8 kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 percent.
Electrical fast transient/burst  IEC 61000-4-4	±2 kV for power supply lines  ±1 kV for input/output lines	±2 kV for power supply lines  ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge  IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV differential mode ±2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions, and voltage variations on power supply input lines  IEC 61000-4-11	<5% $U_T$ (>95% dip in $U_T$ ) for 0.5 cycle  40% $U_T$ (60% dip in $U_T$ ) for 5 cycles  70% $U_T$ (30% dip in $U_T$ ) for 25 cycles  <5% $U_T$ (>95% dip in $U_T$ ) for 5 sec	<5% $U_T$ (>95% dip in $U_T$ ) for 0.5 cycle  40% $U_T$ (60% dip in $U_T$ ) for 5 cycles  70% $U_T$ (30% dip in $U_T$ ) for 25 cycles  <5% $U_T$ (>95% dip in $U_T$ ) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment.  If the user of the CooperSurgical LEEP PRECISION Integrated System requires continued operation during power mains interruptions, it is recommended that the CooperSurgical LEEP PRECISION Integrated System be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field  IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

**NOTE**  $U_T$  is the a.c. mains voltage prior to application of the test level; in this case, 230 V.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environmental Guidance [Notes 1 & 2]
<p>Conducted RF IEC 61000-4-6</p> <p>Radiated RF IEC 61000-4-3</p>	<p>3 Vrms 150 kHz to 80 MHz</p> <p>3 V/m 80 MHz to 2.5 GHz</p>	<p>3 V</p> <p>3 V/m</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the CooperSurgical LEEP PRECISION Integrated System, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = \left[ \frac{3.5}{V_1} \right] \sqrt{P}$ $d = \left[ \frac{3.5}{E_1} \right] \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = \left[ \frac{7}{E_1} \right] \sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ <p>where <math>P</math> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and <math>d</math> is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey<sup>a</sup>, should be less than the compliance level in each frequency range<sup>b</sup>.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

**NOTE 1:** At 80 MHz and 800 MHz, the higher frequency range applies.

**NOTE 2:** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

<sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcasts and TV broadcasts cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the CooperSurgical LEEP PRECISION Integrated System is used exceeds the applicable RF compliance level above, the CooperSurgical LEEP PRECISION Integrated System should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the CooperSurgical LEEP PRECISION Integrated System.

<sup>b</sup> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

## Recommended Separation Distance

### Recommended Separation Distance between portable and mobile RF communications equipment and the CooperSurgical LEEP PRECISION Integrated System.

The CooperSurgical LEEP PRECISION Integrated Systems are intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the CooperSurgical LEEP PRECISION Integrated Systems can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the CooperSurgical LEEP PRECISION Integrated Systems as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter in watts	Separation distance according to frequency of transmitter (Meters) [Notes 1 & 2]		
	150 kHz to 80 MHz $d = \left[ \frac{3.5}{V_1} \right] \sqrt{P}$	80 MHz to 800 MHz $d = \left[ \frac{3.5}{E_1} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[ \frac{7}{E_1} \right] \sqrt{P}$
0.01	0.1167	0.1167	0.2333
0.1	0.3689	0.3689	0.7379
1	1.1667	1.1667	2.3333
10	3.6894	3.6894	7.3789
100	11.667	11.667	23.333

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output rating of the transmitter in watts (W) according to the transmitter manufacturer.

For the CooperSurgical LEEP PRECISION Integrated System:

$$V_1 = 3 \text{ Vrms}$$

$$E_1 = 3 \text{ V/m}$$

**NOTE 1:** At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

**NOTE 2:** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

## Section 17 Liability Statement

CooperSurgical guarantees the safety, reliability, and performance of the LEEP PRECISION Integrated System only if the installation, recalibrations, and repairs are performed by personnel authorized by CooperSurgical and if it is used in compliance with given instructions in an area that meets all applicable IEC requirements.

## Section 18 Warranty

CooperSurgical, Inc., warrants that the LEEP PRECISION Integrated System (the "Product") will be free from defects in materials and workmanship for a period of one (1) year from the original date of purchase.

If the Product should become inoperable due to a defect in material or workmanship during this one year warranty period, CooperSurgical will, at its option, repair or replace the Product. This limited warranty does not include replacement or service to repair damage resulting from improper installation, external electrical fault,

accident, disaster, use for a purpose other than that for which originally designed or indicated in this manual, negligence, modification, service or repair by personnel not authorized by CooperSurgical, or normal wear and tear, and also does not apply to disposable or single- or limited-use items or components. The sole and exclusive remedy under this limited warranty shall be repair or replacement as provided herein.

The foregoing limited warranty states the sole warranty made by CooperSurgical with respect to the Product and all parts thereof, and is in lieu of any other warranty by CooperSurgical with respect to Product.

COOPERSURGICAL NEITHER MAKES NOR GRANTS ANY OTHER WARRANTY, EITHER EXPRESS OR IMPLIED, WITH RESPECT TO THE PRODUCT, INCLUDING WITHOUT LIMITATION, WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL COOPERSURGICAL BE LIABLE FOR ANY DAMAGES ARISING OUT OF THE LOSS OF USE OF THE PRODUCT, OR ANY OTHER INCIDENTAL OR CONSEQUENTIAL DAMAGES, WHETHER OR NOT COOPERSURGICAL HAS ADVANCE KNOWLEDGE OF THE POSSIBILITY OF SAME.

No person, agent, distributor, dealer, or company is authorized to change or modify the terms of this Limited Warranty.

**Only CooperSurgical is authorized to service or repair this unit. Do not disassemble the device. There are no user-serviceable components within the housing.**

## Section 19 Service and Repair

**LEEP PRECISION Generator:** There are no customer-serviceable parts.

**LEEP PRECISION Cart:** There are no customer-serviceable parts.

**LEEP PRECISION Integration Unit:** There are no customer-serviceable parts.

**LEEP PRECISION Smoke Evacuator:** There are no customer-serviceable parts (excluding intermittent filter changes).

Only CooperSurgical, Inc., is authorized to service or repair the LEEP PRECISION Integrated System. If repair is attempted outside the factory, the warranty is considered void. CooperSurgical is not responsible for any injury resulting from repairs made by other individuals or organizations not certified by CooperSurgical, Inc. If a repair is needed, the equipment *must be sanitized* before it is returned to the factory and carefully packaged in a protective carton.

**On the note inserted in the box, please provide the following information:**

- Customer and contact information on the Repair Authorization Form (download from CooperSurgical website) or on company letterhead
- Nature of the problem
- Description of the item returned
- Serial number (if applicable)

All shipments must be prepaid. COD packages will not be accepted. Return carton to:

### **CooperSurgical, Inc.**

Attention: Repair Department  
95 Corporate Drive • Trumbull, CT 06611 USA  
Phone: +1 (203) 601-5202 • (800) 444-8456  
Fax: +1 (203) 601-4743

# Section 20 Explanation of Symbols

-  Reorder number
-  Serial number
-  Consult instructions for use
-  Caution
-  Non-sterile
-  **APPLIED PART SYMBOL** = Electrode that is connected to a patient's body.
-  **DEFIBRILLATOR PROOF SYMBOL** = indicates that the device will not be damaged if the defibrillator is active.
- BF-TYPE EQUIPMENT** = Type of equipment is a B-Type equipment with an F-Type applied part.  
**B-Type** equipment is a piece of equipment providing a particular degree of protection against electric shock, particularly regarding:
  - Allowable leakage current
  - Reliability; of the protective earth connection.**F-Type** = Isolated applied part. Applied part isolated from all other parts of the equipment to such a degree that the patient leakage current allowable in single fault condition is not exceeded when a voltage equals to 1.1 times the highest rated mains voltage is applied between the applied part and earth.
-  Dangerous voltage
-  Fuse, Replace only with Specified Type and Ratings
-  Alternating Current
-  Equipotentiality Symbol (wires carry the same voltage)

-  Earth (Ground)
-  Medical equipment with respect to electrical shock, fire and mechanical hazard only in accordance with UL60601-1 and CAN/CSA C22.2 No.601.1
-  Date of Manufacture
-  Manufacturer
- Rx Only**  
**Caution:** U.S. Federal law restricts this device to sale by or on the order of a physician.
-  This side up
-  Keep Dry
-  Fragile
-  In order to preserve, protect and improve the quality of the environment, protect human health and utilize natural resources prudently and rationally – do not dispose of waste electrical or electronic equipment (WEEE) as unsorted municipal waste. Contact local WEEE disposal sites.

Velcro® is a registered trademark of Velcro Industries B.V.  
 © 2024 CooperSurgical, Inc.  
 Made in the USA

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Teléfono: +1 (203) 601-9818  
Fax: +1 (203) 601-4747

CooperSurgical, Inc.  
95 Corporate Drive  
Trumbull, CT 06611 EE. UU.

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VELCRO® es una marca registrada de Velcro Industries B.V.