# Ring Pessary PVC

wallace®

en

INSTRUCTIONS FOR USE

fr

MODE D'EMPLOI

de

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it

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es

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pt

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nl

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sv

BRUKSANVISNING

da

BRUGSANVISNING

no

BRUKSANVISNING

fi

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el

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EC REP

EMERGO EUROPE  
Prinsessegracht 20  
2514 AP The Hague  
The Netherlands

CooperSurgical, Inc.  
95 Corporate Drive  
Trumbull, CT 06611 USA

These instructions are for use with the following Wallace® products:

Catalog No. 700/300/050 - 700/300/110 *Ring Pessary PVC*

These instructions contain important information for safe use of the product. Read the entire contents of these Instructions for Use, including Warnings and Cautions, before using this product. Failure to properly follow warnings, cautions and instructions could result in death or serious injury to the patient and/or clinician.

NOTE: DISTRIBUTE THIS INSERT SHEET TO ALL PRODUCT LOCATIONS.

#### DESCRIPTION

- A range of Ring Pessaries manufactured in PVC for the palliative treatment of uterine displacement and stress incontinence.
- The PVC Ring Pessary is flexible with a large cross section.
- Wallace® Ring Pessaries are available on the UK National Health Prescription.

#### WARNINGS

- Insertion and removal of Ring Pessaries should only be performed by competent trained personnel.
- After removal of the Ring Pessary an inspection of the vagina should be performed to check for ulceration, bleeding, infection.
- This product contains DEHP – data derived from animal and epidemiological studies outline the following risks, but are not limited to: Adverse effects on fetal development, kidney and liver toxicity, negative impact on fertility. Special attention on the application of phthalate products should be given for the following patient populations: pediatrics, pregnant and nursing women.

#### INDICATIONS FOR USE

- Wallace Ring Pessaries provide vaginal support to prevent the uterus from descending.
- Designed to ease uterine prolapse.
- Ring Pessaries are also used to correct other uterine displacements and to alleviate stress incontinence by pressure through the anterior vaginal wall onto the urethra.

#### CONTRAINDICATIONS

- Ring Pessaries are not suitable for patients with active infections of the pelvis or vagina.
- Ring Pessaries should not be prescribed for patients who are noncompliant, or unlikely to follow-up.

#### PRECAUTIONS

- Inspection of the vagina using a speculum is recommended prior to insertion or replacement of a Ring Pessary.
- The Wallace® Ring Pessaries are supplied non-sterile and are reusable for single patient use only. Products must be disposed of after 6 months.

#### INSTRUCTIONS FOR USE

- 1 Determine the size of Ring Pessary required by currently accepted medical techniques.
- 2 The Ring Pessary may be immersed in hot water to warm it and increase its flexibility if required.
- 3 Wearing sterile gloves, hook the middle finger around the outer rim of the Ring Pessary in the center until the edges make contact.
- 4 To ease insertion, a water-soluble lubricating jelly may be applied to the front edge of the Ring Pessary. Holding the Ring Pessary vertically, insert into the vagina angling it down towards the floor of the pelvic basin.
- 5 Remove fingers and thumb and push the Ring Pessary into place ensuring it is on the bony ledge behind the symphysis pubis.
- 6 Check the position of the Ring Pessary, the cervix should be felt protruding through it. The Ring Pessary should be across the vagina with the back edge in the posterior fornix and the front edge behind the symphysis pubis.

#### INSTRUCT/ADVISE PATIENT

- 1 Instruct the patient that the Ring Pessary should only be inserted and removed as directed by the attending physician.
- 2 Instruct the patient that, after removal, she should wash the Ring Pessary with soap and water and rinse before re-insertion.
- 3 Instruct the patient that it is important for her to attend follow-up appointments.
- 4 Advise the patient that there is a possibility that she may experience discomfort or more vaginal discharge than usual and instruct her to report any such side effects to her physician immediately.
- 5 Advise the patient how long she should be expected to use the Ring Pessary.
- 6 Advise the patient that the Ring Pessary may interfere with sexual intercourse. If any discomfort is experienced, her physician should be consulted.

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Wallace® is a registered trademark of CooperSurgical, Inc. and its subsidiaries. CooperSurgical is a registered trademark of CooperSurgical, Inc.