

Anospec® | by CooperSurgical®

# A BRILLIANT SOLUTION



Cordless high-powered illumination for 30+ minutes that doesn't generate damaging heat.



Save time and increase efficiency – ready to use with integrated LED light source and obturator.



Single-use speculum reduces cross contamination and eliminates reprocessing costs.




The device's clear material allows for 360 degree visibility.

## Ordering Information

| NUMBER  | DESCRIPTION                       |
|---------|-----------------------------------|
| C060100 | Anospec® (90 mm x 18 mm Slotted)  |
| C060110 | Anospec® (96 mm x 23 mm Slotted)  |
| C060120 | Anospec® (103 mm x 18 mm Beveled) |

Batteries included in each unit. Sold clean, non-sterile.



[Learn More](#)

**IMPORTANT SAFETY INFORMATION ANOSPEC® Single-Use** Anoscope is a lighted anal speculum intended for physician use to examine the anal sphincter, anus, rectum, and, using additional devices, to perform various diagnostic and therapeutic procedures. **Warnings: Never** use with laser equipment, gamma radiation, autoclave sterilization. Refer to Instructions for Use, for detailed instructions and risk information. [www.coopersurgical.com/anospec-instructions-for-use](http://www.coopersurgical.com/anospec-instructions-for-use)

800.243.2974 | [www.coopersurgical.com](http://www.coopersurgical.com)



CooperSurgical®