

Product Specifications

LEEP PRECISION Integrated System

Dimensions: (W x D x H)
22.44 in. x 22.07 in. x 41.22 in.
57.0 cm x 56.06 cm x 104.70 cm

Weight: 88 pounds, 11 ounces (40.23 kg)

Input Voltage: 120 VAC, 60Hz

VA Rating: 132 VA

Fuses: 10A, 250V, T-Type, Slow Blow

LEEP PRECISION Generator

Dimensions: (W x D x H)
9.38 in. x 11.63 in. x 4.25 in.
23.8 cm x 29.5 cm x 10.8 cm

Weight: 10 pounds, 9 ounces (5.25 kg)

Input Voltage: 120 VAC, 60Hz

Maximum Current: 3.15 amps

Power Line Leakage: Less than 50 microamps

Fuses: Two T 3.15A (5 x 20 mm) 250V

LEEP PRECISION Smoke Evacuator

Dimensions: (W x D x H)
8½ in. x 9 in. x 22 in.
216 cm x 229 cm x 559 cm

Weight: 24 pounds (11 kg)

Air Flow: 35 CFM minimum

Filter: Pre-Filter, charcoal filter for odor removal and hydrophobic ULPA filter rated at 99.999% efficiency level for 0.014 micron particles*.

Input Voltage: 120 VAC, 60 Hz

Fuses: 15 amps, 250V, T-Type

LEEP PRECISION Cart

Dimensions: (W x D x H)
22.44 in. x 22.07 in. x 41.22 in.
57.0 cm x 56.06 cm x 104.70 cm

Weight: 54 pounds, 2 ounces (24.55 kg)

LEEP PRECISION™ Integrated System Ordering Information

Code	Product
LEEP PRECISION Integrated System	
LP-10-120	LEEP PRECISION Integrated System, 120V
LEEP PRECISION Generator	
LP-20-120	LEEP PRECISION Generator, 120V
LEEP PRECISION Smoke Evacuator	
LP-30-120	LEEP PRECISION Smoke Evacuator, 120V, stand-alone
LEEP PRECISION Cart	
LP-40-001	LEEP PRECISION Cart
LEEP PRECISION Generator Disposables	
LP-50-101	Hand Switch-operated Push-Button Cautery Pencils (10/Box)
LP-50-102	Foot Pedal-Operated Cautery Pencils (10/Box)
LP-50-103	Foot Pedal Adapter (1/Bag)
LP-50-201	Dual Dispersive Electrodes with Connector (10/Box)
LP-50-202	Dual Dispersive Electrodes without Connector (connect to Reusable Cable) (10/Box)
LP-50-203	Reusable Cable (1/Bag)
LEEP PRECISION Smoke Evacuator Disposables	
52560	Smoke Evacuator Pre-Filter, Single-Patient Use (10/Bag)
6082	Smoke Evacuator ULPA Filter (1/Box)
6083	Smoke Evacuator Reducer, Single-Patient Use (10/Bag)
6084	Smoke Evacuator Tubing – 6 ft. long, ½ in. OD (Sterile, Single Patient Use) (10/Box)



CooperSurgical

CooperSurgical is the leading company dedicated to providing medical devices and procedure solutions that improve health care delivery to women regardless of clinical setting. Our company is fostering that position through expansion of its core businesses and introduction of advanced technology-based products, which aid clinicians in the management and treatment of commonly seen conditions.

800.243.2974 | 203.601.5200 | coopersurgical.com | 95 Corporate Drive | Trumbull, CT 06611
Also on [YouTube.com/coopersurgical](https://www.youtube.com/coopersurgical)

LEEP PRECISION™ Integrated System

Office Diagnostics and Procedures

The Preferred Choice for Functionality, Performance, and Precision for Loop Electrosurgical Excision Procedures

- Hand or foot activation for operating precision
- Smoke evacuator eliminates odor at 99.999 percent efficiency level
- Fully integrated LEEP station with storage
- Safe and easy to operate
- Heavy-duty casters facilitate mobility
- Elegantly designed for the modern procedure room



CooperSurgical

The Preferred Choice for Functionality, Performance and Precision for LEEP Procedures

LEEP PRECISION™ Integrated System has been designed to meet the unique requirements of the clinic and practice-based facility. Fully integrated and mobile, the LEEP PRECISION Integrated System combines three components – LEEP PRECISION Cart, LEEP PRECISION Generator, and the LEEP PRECISION Smoke Evacuator – for maximum efficiency



LEEP PRECISION Cart



- Generator and smoke evacuator in one compact station for added mobility and functionality
- Designed to work in small exam rooms
- Heavy-duty casters ensure easy mobility, two locks to secure placement
- Two convenient storage shelves
- Comfort-grip handle for easy movement and cart positioning
- Cable hooks to manage cable wires before, during, and after procedure
- Elegant design for modern medical facilities

LEEP PRECISION Generator



- LED displays located on the front for precise power selection, delivery, and ease of use
- Flush faceplate membrane facilitates operation and cleaning
- Microprocessor controlled for increased precision, accuracy, reproducibility, and safety
- Provides a choice of Cut, Blend, and Coag wave forms to accommodate subtle differences in technique and electrode performance
- Hand or foot activation options for easy operating precision
- Patient return monitor disables power automatically if grounding pad becomes loose
- Self-monitors all detectors



LEEP PRECISION Smoke Evacuator

- High air flow efficiently captures smoke plumes (35 CFM min.)
- Triple-stage filtration captures airborne particulate matter, vapor, and odor with a 99.999 percent efficiency level
- Virtually maintenance-free
- Replacement filters available