



Premier Products for Clinicians Treating Mom and Baby

Labor and Delivery Products Catalog







Products for obstetrical and neonatal care



Physicians have long depended on our instruments throughout the labor and delivery process.

Our portfolio for obstetrical and neonatal care includes devices for vacuum-assisted delivery, fetal monitoring, amniocentesis, PROM testing, infant critical care and C-section.

Labor and Delivery Brands

At CooperSurgical, we develop and market innovative solutions for labor and delivery procedures. Within this catalog, you will find an extensive range of instruments and devices that can be used in your practice to improve patient care.



Fetal Pillow®

INSORB® Subcuticular Skin Stapler

Actim® PROM

MityOne®

Mystic[®] II

Mityvac[®]

GUARDIAN Vaginal Retractor®

MityHook™

INCA®

TransWarmer[®]

WarmGel®

C SAFE®

OB/Mobius®

Table of Contents


Fetal Pillow [®]
Subcuticular Skin Stapler
POC Test for Ruptured Membranes Actim® PROM
Vacuum-Assisted Delivery Systems MityOne®
Mystic® II and Mityvac®
Comfort Belly Bands
Episiotomy Retractors
GUARDIAN Vaginal Retractor®

Neonatal CareINCA® Infant Nasal CPAP Assembly.14TransWarmer® Infant Transport Mattress.15WarmGel® Infant Heel Warmer.15
C-Section Devices C SAFE®
Warranty Information/ Customer Service Information

CE Available for sale where CE mark is required

CA* Available for sale in Canada





^{*} Please refer to IFU for detailed information



Actim® PROM

Accurate PROM Detection, Even in the Presence of More Than Trace Blood and Other Contaminants¹



Actim® PROM

NUMBER	DESCRIPTION
30831ETUS	Actim PROM 10 Test Kit
30800ETUS	Actim PROM Control Test Kit

1. The test has been designed to minimize interference from bleeding, but in cases of heavy bleeding the blood locally may have a higher concentration of IGFBP-1 protein. In these cases, a positive result should be interpreted with caution. Actim PROM IFU.



MityOne®


Vacuum Pumps Self-limiting M-Style® vacuum Mushroom® Cup · Shorter stem Freely rotating stem Prominent vacuum release positioned closer to operator Ergonomic design Semi-rigid MitySoft® Bell Cup MityOne® Vacuum Pumps CA ★ CE NUMBER **DESCRIPTION** 10067 MityOne Pump with M-Style® Mushroom® Cup (12/Box, Min. 12) 10068 MityOne Pump with MitySoft® Bell Cup (12/Box, Min. 12)

Mystic II and Mityvac

Vacuum-Assisted Delivery Pumps and Accessories


One-Piece Mystic II System





The Original Mityvac® Two-Piece System

The original two-piece Mityvac system is comprised of the reusable Mityvac pump and a variety of cup options.



Mityvac° and Mystic° II Vacuum-Assisted Delivery Pumps and Accessories

16177 Mityvac Reusable Pump Replacement Gauge (1/Box) 16710 Mityvac Reusable Pump Maintenance and Repair Kit (1/Box) 16178 Mityvac Disposable Tubing and Filter – Use with Mityvac Reusable Cup (12/Box) 10502 MitySoft Reusable Cup Replacement Valve (1/Box) 10057 Mystic II Pump with M-Style® Mushroom® Cup (12/Box) 10058 Mystic II Pump with MitySoft® CA*			
16177 Mityvac Reusable Pump Replacement Gauge (1/Box) 16710 Mityvac Reusable Pump Maintenance and Repair Kit (1/Box) 16178 Mityvac Disposable Tubing and Filter – Use with Mityvac Reusable Cup (12/Box) 10502 MitySoft Reusable Cup Replacement Valve (1/Box) 10057 Mystic II Pump with M-Style® Mushroom® Cup (12/Box) 10058 Mystic II Pump with MitySoft® CA*	NUMBER	DESCRIPTION	
(1/Box) 16710 Mityvac Reusable Pump Maintenance and Repair Kit (1/Box) 16178 Mityvac Disposable Tubing and Filter – Use with Mityvac Reusable Cup (12/Box) 10502 MitySoft Reusable Cup Replacement Valve (1/Box) 10057 Mystic II Pump with M-Style® CA* Mushroom® Cup (12/Box) 10058 Mystic II Pump with MitySoft® CA*	10022	Mityvac Reusable Handheld Pump (1/Box)	CA* (€
Repair Kit (1/Box) 16178 Mityvac Disposable Tubing and Filter – Use with Mityvac Reusable Cup (12/Box) 10502 MitySoft Reusable Cup Replacement Valve (1/Box) 10057 Mystic II Pump with M-Style® Mushroom® Cup (12/Box) 10058 Mystic II Pump with MitySoft® CA*	16177	, , , , ,	
Use with Mityvac Reusable Cup (12/Box) 10502 MitySoft Reusable Cup Replacement Valve (1/Box) 10057 Mystic II Pump with M-Style® CA® Mushroom® Cup (12/Box) 10058 Mystic II Pump with MitySoft® CA®	16710		
(1/Box) 10057 Mystic II Pump with M-Style® CA* Mushroom® Cup (12/Box) 10058 Mystic II Pump with MitySoft® CA*	16178	, ,	
Mushroom® Cup (12/Box) 10058 Mystic II Pump with MitySoft® CA*	10502	1 1	
, ,	10057	, ,	CA * (€
Бен Сар (12/Вох)	10058	Mystic II Pump with MitySoft® Bell Cup (12/Box)	CA * (€

Mityvac Vacuum-Assisted Delivery Cups ca∗ <€

NUMBER	DESCRIPTION
10004	Standard Pearl Edge™ Extraction Cup (12/Box)
10007LP	M-Style® Mushroom® Cup (50 mm) with 4' Tubing and Filter (12/Box)
10008	Super M-Style Mushroom Cup (56 mm) with 4' Tubing and Filter (12/Box)
10020	MitySoft® Bell Cup (64 mm) with 4' Tubing and Filter (12/Box)
10137	M-Select® Mushroom Cup (50 mm) with 4' Tubing and Filter (12/Box)
10500	Reusable Vacuum Delivery Cup (70 mm) (1/Box)

Comfort Belly Bands

Abdominal Fetal Monitoring Band

Belly Bands are available in a wide range of widths.

- Breathable material
- Flexible precut lengths available in 30-foot rolls or individually wrapped packages of precut 14-inch lengths



Comfort Belly Bands

NUMBER	DESCRIPTION
13010	10" Width Belly Bands on Roll (10" x 10 yds) (1 Roll)
13011	10" Width Belly Bands – Precut (10" x 14") (25/Box)
13012	12" Width Belly Bands on Roll (12" x 10 yds) (1 Roll)
13013	12" Width Belly Bands – Precut (12" x 14") (25/Box)
13014	14" Width Belly Bands on Roll (14" x 10 yds) (1 Roll)
13015	14" Width Belly Bands – Precut (14" x 14") (25/Box)
13016	16" Width Belly Bands on Roll (16" x 10 yds) (1 Roll)
13017	16" Width Belly Bands – Precut (16" x 14") (25/Box)
13019	18" Width Belly Bands – Precut (18" x 14") (25/Box)





GUARDIAN Vaginal Retractor

Episiotomy Retractor

The GUARDIAN Vaginal Retractor provides visualization and access to the complete circumference of the vagina and perineal tissue during episiotomy or laceration repair.

NUMBER	DESCRIPTION
64-330	Episiotomy Retractor (6/Box)



$MityHook^{\mathsf{m}}$

Amniotome - MityHook™ CA*

For amniotic membrane rupture.

NUMBER	DESCRIPTION
10035	MityHook, Amniotome, Sterile, 9" (100/Box)





INCA® Infant Nasal CPAP Assembly

Neonatal Nasal CPAP



INCA® Complete Set CA ★ CE

Consists of a nasal cannula, two lengths of Stay-Flex tubing, two tubing adapters, one in-line "T" pressure adapter, a knit cap with hook/loop strap, two foam blocks and a hook/loop chin strap. One INCA Sizing Gauge is included in every box of five sets.

NUMBER	DESCRIPTION
44-2707	Set with 7.5 Fr Cannula (5/Box)
44-2709	Set with 9 Fr Cannula (5/Box)
44-2710	Set with 10.5 Fr Cannula (5/Box)
44-2712	Set with 12 Fr Cannula (5/Box)
44-2715	Set with 15 Fr Cannula (5/Box)



INCA Accessory Pack 🧯

Each pack includes five knit caps with hook/loop strap, five chin straps with ten hook fasteners (not attached) and ten foam blocks.

NUMBER	DESCRIPTION
44-2001	X-Small Pack (5 XS Caps, 10 Foam Blocks)
44-2002	Small Pack (5 S Caps, 10 Foam Blocks)
44-2003	Medium Pack (5 M Caps, 10 Foam Blocks)
44-2004	Large Pack (5 LG Caps, 10 Foam Blocks)



INCA Replacement Set CA ★ CE

Includes the nasal cannula, two lengths of Stay-Flex tubing, two tubing adapters and one in-line pressure "T" adapter.

NUMBER	DESCRIPTION
44-0707	Set with 7.5 Fr Cannula (5/Box)
44-0709	Set with 9 Fr Cannula (5/Box)
44-0710	Set with 10.5 Fr Cannula (5/Box)
44-0712	Set with 12 Fr Cannula (5/Box)
44-0715	Set with 15 Fr Cannula (5/Box)



INCA Nasal Sizing Gauge ← €

For appropriate nasal prong fitting.

NUMBER	DESCRIPTION
025-001	INCA Nasal Sizing Gauge (1/Pack)



NEO-prep[™] Skin Barrier <€

Neo-prep Skin Barriers are precut and non sterile.

NUMBER	DESCRIPTION
44-8000	Precut Hypoallergenic Skin Barrier, non sterile, small (10/Box)
44-8001	Precut Hypoallergenic Skin Barrier, non sterile, large (10/Box)



NEO-fit™ Neonatal Endotracheal Tube Grip

Facilitates Endotracheal Tube Management The NEO-fit provides endotracheal tube stabilization. Designed to firmly grip a neonatal endotracheal tube, NEO-fit utilizes metal clips to mechanically grip the soft material of the tube.

NUMBER	DESCRIPTION
42-2540	NEO-fit Neonatal Endotracheal Tube Grip Set (20/Box). Each Set includes: (1) plastic tube grip, with hook/loop strap, (2) foam foot pads, (2) foam "lollipop" tapes




TransWarmer[®] Infant Transport Mattress

TransWarmer® Infant Transport Mattress

NUMBER DESCRIPTION

20421 TransWarmer Infant Transport Mattress (6/Box)

 \in







WarmGel® Infant Heel Warmer

WarmGel® Infant Heel Warmer

Gel-filled, disposable heel warmer.

NUMBER	DESCRIPTION	
20418	WarmGel Infant Heel Warmer (100/Box)	C€





C-Section Devices

C SAFE®

Uterine Safety Scalpel

Designed to help reduce the risk of fetal lacerations during cesarean section procedures.

C SAFE® Uterine Safety Scalpel

NUMBER	DESCRIPTION
BX101	C SAFE Safety Scalpel (10/Box)






C-Section Devices


OB/Mobius®

Self-Retaining Elastic Retractor

Available in two sizes for C-sections.

Smaller diameter Mobius retractors are also available for other laparotomy and minilaparotomy procedures including hysterectomy, cystectomy, oophorectomy and more.





OB/Mobius® CA₩ (€

NUMBER	MBER DESCRIPTION	
900-517	Medium-Large Retractor, 17 cm l.D. ring, for incisions to 15 cm (5/Box)	
900-520	Large Retractor, 22 cm I.D. ring, for incisions to 17 cm (5/Box)	

Also available in smaller sizes for minilaparotomies.

OR/Mobius® CA* CE

NUMBER	DESCRIPTION
900-508	Mini Retractor, 8 cm I.D. ring, for incisions to 6 cm (5/Box)
900-500	Small Retractor, 11 cm l.D. ring, for incisions to 9 cm (5/Box)
900-515	Medium Retractor, 15 cm I.D. ring, for incisions to 14 cm (5/Box)



Committed to Our Customers

At CooperSurgical, we strive for the highest level of expertise and knowledge of our products and to develop an ongoing relationship with each of our customers that is built on trust, quality and consistency. We look to provide service that exceeds not only our customers' expectations, but industry standards as well. CooperSurgical's Customer Service Department provides full–service order entry with extended hours to facilitate order placement. As an added convenience we also now offer online order entry. Call now to have your Customer Service Representative set up your account today. Once you are registered for online ordering, you will be able to place orders, track shipments, review prior purchase history and schedule orders for future delivery.

CE Available for sale where CE mark is required

CA* Available for sale in Canada

CE and CA labeling is subject to change. Contact Customer Service or your local sales representative for the most up-to-date information.

5 Ways to Order



Local Sales Representative

www.coopersurgical.com/pages/salesrepsearch.aspx



Email

U.S.: Orders@coopersurgical.com

Int'l: InternationalCustomerservice@coopersurgical.com



E-Ordering (U.S. Only)

Call 800.243.2974 to set up your account today.



Phone Orders

800.243.2974 (8:00 AM - 8:30 PM EST)

Int'l: 203.601.9818



Fax Orders

800.262.0105 (24 hours a day)

Int'l: 203.601.4747

Warranty Information

CooperSurgical, Inc. ("CooperSurgical") warrants that its products offered in this catalog, **other than** single-use products and consumables, will be free from defects in materials and workmanship for one year from the date of purchase. CooperSurgical warrants that its single-use products and consumables offered in this catalog will, through the earlier of 90 days from the date of purchase of each such product and the expiration date of that product stated on that product's package, be free from defects in materials and workmanship that are discovered prior to or during the single use of that product.

If CooperSurgical determines that a product offered in this catalog fails to conform to the CooperSurgical warranty applicable to that product, CooperSurgical will, as the sole remedy for that failure to conform, repair or replace that product, at CooperSurgical's discretion, free of charge.

To return a product to CooperSurgical, a customer must comply with CooperSurgical's Returned Goods Policy described in this catalog. A customer will not have any remedy if a CooperSurgical product does not conform to the warranty for that product unless that product is returned to CooperSurgical in accordance with that Returned Goods Policy. CooperSurgical will ship returned products that it repaired or replaced under warranty to the customer who returned those products, at CooperSurgical's expense F.O.B. the customer's facility. Under all other circumstances, CooperSurgical will ship returned products to the customer who returned those products at the customer's expense F.O.B. CooperSurgical's facility.

CooperSurgical's warranties, including those in this catalog, do not cover damage caused by misuse, improper care, improper use of chemicals or cleaning methods, loss, theft, servicing or sharpening by non-authorized personnel or negligent or intentional conduct on the part of the owner or user of the product. CooperSurgical's warranties do not cover normal wear and tear or general maintenance, including cleaning and sharpening services required as a result of normal use. Any modifications or changes to a product will void that product's warranty.

THE ABOVE WARRANTIES ARE IN LIEU OF, AND COOPERSURGICAL HEREBY DISCLAIMS, ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, WRITTEN OR ORAL, WITH RESPECT TO COOPERSURGICAL'S PRODUCTS, INCLUDING THE WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. NO TERMS, CONDITIONS, UNDERSTANDINGS OR AGREEMENTS THAT PURPORT TO MODIFY THE ABOVE WARRANTIES OR THAT MAKE ANY ADDITIONAL WARRANTIES FOR ANY COOPERSURGICAL PRODUCT SHALL HAVE ANY LEGAL EFFECT UNLESS MADE IN WRITING AND SIGNED BY AN AUTHORIZED COOPERSURGICAL CORPORATE OFFICER.

COOPERSURGICAL SHALL NOT UNDER ANY CIRCUMSTANCES BE LIABLE FOR LOST PROFITS, DAMAGES FROM LOSS OF USE OR LOST DATA, OR INDIRECT, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES UNDER ITS WARRANTIES OR OTHERWISE FOR ANY CLAIM RELATED TO COOPERSURGICAL'S PRODUCTS, EVEN IF COOPERSURGICAL HAS BEEN ADVISED, KNEW OR SHOULD HAVE KNOWN OF THE POSSIBILITY OF SUCH DAMAGES. COOPERSURGICAL'S LIABILITY WITH RESPECT TO A PRODUCT COVERED BY A WARRANTY OR OTHERWISE SHALL BE LIMITED IN ALL CIRCUMSTANCES TO THE PURCHASE PRICE OF THAT PRODUCT.

Returned Goods Policy

- 1. Goods will be accepted for return for the following reasons:
 - If shipment was made without Customer's authorization or order
 - If incorrect goods were shipped
 - If defective goods are covered by the standard warranty
- 2. To return goods, you must contact a Customer Service Representative by telephone (1.800.243.2974 or 1.203.601.5200) for a Returned Merchandise Authorization (RMA) number. Items will not be accepted without a RMA number.

Please have the following information:

- Reason you wish to return the goods
- · Quantity, description, part number, serial number of the goods
- · Date of Receipt of Order
- Customer's Purchase Order and the CooperSurgical Invoice Number All used instruments must be cleaned and sterilized prior to being shipped for return.
- 3. Shipment must be sent prepaid by Customer. Freight collect shipments will not be accepted, and goods will be returned to sender.
- **4.** If Customer intends to return equipment ordered in error, the following restocking charges and terms will apply:
 - There is no restocking fee if goods are returned within 30 days of shipment
 - There is a restocking fee of 25% if goods are returned within 31-60 days of shipment
 - Goods must be returned unused, in the original carton, and in marketable condition
 - Refurbishing and replacement charges will be added to the restocking charges for damaged or missing items
 - No return after 60 days
 - No refund on sterile, single-use disposable products

Send All Returns to: CooperSurgical

95 Corporate Drive Trumbull, CT 06611

Trademark Notices

The marks bearing the symbols "TM" and "®" are trademarks and registered trademarks of CooperSurgical, Inc.

Actim® PROM is a registered trademark of Medix Biochemica.



Fetal Pillow®

INSORB® Subcuticular Skin Stapler

Actim® PROM

MityOne®

Mystic[®] II

Mityvac[®]

GUARDIAN Vaginal Retractor®

MityHook™

INCA®

TransWarmer®

WarmGel®

C SAFE®

OB/Mobius®





800.243.2974 www.coopersurgical.com 95 Corporate Drive, Trumbull, CT 06611